



# eBABYHR Mounting Instructions for 553KL/553KG



PELSTAR, LLC 9500 West 55<sup>th</sup> St. McCook, IL 60525-7110 USA

[www.homscales.com](http://www.homscales.com)

© Pelstar, LLC 2018

weigh**easier**®






# Health o meter®

Professional

Thank you for purchasing this Health o meter® Professional product. Please read these instructions carefully and keep it for easy reference or training.

## Notes:

- Ensure that the scale is off during the height rod mounting process. If the scale was powered on during the process, be sure to turn the scale off and power on again before weighing the first patient.
- If the height rod needs to be removed, ensure the scale is off before removal and off before replacing the height rod onto the scale.

<p>1. Gently place the scale upside down on a flat, stable surface. Using a soft cloth and Isopropyl alcohol, clean the surface around the handhold where the bracket will be mounted, as shown.</p>	
<p>2. Obtain the mounting bracket and remove the backing of the double sided tape. Align the bracket onto the scale around the handhold, as shown. Push down to affix the tape and bracket onto the scale. <b>Warning:</b> The tape is very strong; do not push down until bracket is properly positioned.</p>	
<p>3. Obtain the four self-drilling screws from the packaging. Using a power drill, insert the four screws through the holes in the bracket into the underside of the scale.</p>	
<p>4. Tip the scale onto its side. Obtain the eBABYHR. Align the mounting channel slots on the back of the height rod display with the bracket. Push down on the oval tab and slide the height rod's mounting channel onto the bracket until you hear a click.</p>	
<p>5. Obtain the two thumb screws from the packaging. Insert the thumb screws through the holes in the bracket and into the back of the height rod display. Hand tighten the screws. Ensure that the stationary paddle of the height rod is flush with the left side of the weighing tray. Mounting is now complete.</p>	

# WARRANTY

## **Limited Warranty**

### **What does the Warranty Cover?**

This Health o meter<sup>®</sup> Professional product is warranted from date of purchase against defects of materials or in workmanship for a period of two (2) years. If product fails to function properly, return the product, freight prepaid and properly packed to Pelstar, LLC (see "To Get Warranty Service", below, for instructions). If the manufacturer determines that a defect of material or in workmanship exists, the customer's sole remedy will be replacement of the product at no charge. Replacement will be made with a new or remanufactured product or component. If the product is no longer available, replacement may be made with a similar product of equal or greater value. All replaced parts are covered only for the original warranty period.

### **Who is Covered?**

The original purchaser of the product must have proof of purchase to receive warranty service. Please save your invoice or receipt. Pelstar dealers or retail stores selling Pelstar products do not have the right to alter, or modify or any way change the terms and conditions of this warranty.

### **What is Excluded?**

**Your warranty does not cover normal wear of parts or damage resulting from any of the following:** negligent use or misuse of the product, use on improper voltage or current, use contrary to the operating instructions, abuse including tampering, damage in transit, or unauthorized repair or alternations. Further, the warranty does not cover natural disasters, such as fire, flood, hurricanes and tornadoes. This warranty gives you specific legal rights, and you may also have other rights that vary from country to country, state to state, province to province or jurisdiction to jurisdiction.

**To get Warranty Service make sure you keep your sales receipt or document showing proof of purchase.** Call (+1) 800-638-3722 or (+1) 708-377-0600 to receive a return authorization (RA) number, which must be included on the return label. Attach your proof of purchase to your defective product along with your name, address, daytime telephone number and description of the problem. Carefully package the product and send with shipping and insurance prepaid to:

Pelstar, LLC  
Attention R/A# \_\_\_\_\_  
Return Department  
9500 West 55<sup>th</sup> Street  
McCook, IL 60525



**PELSTAR, LLC**  
**9500 West 55th St. McCook, IL 60525-7110 USA**  
**1-800-638-3722 or 1-708-377-0600**

PLEASE REGISTER YOUR PRODUCT FOR WARRANTY COVERAGE AT:  
[www.homscales.com](http://www.homscales.com)

Health o meter<sup>®</sup> is a registered trademark of Sunbeam Products, Inc. used under license. Health o meter<sup>®</sup> Professional products are manufactured, designed, and owned by Pelstar, LLC. We reserve the right to improve, enhance, or modify Health o meter<sup>®</sup> Professional product features or specifications without notice.

© Pelstar, LLC 2018