


Health o meter[®]
P r o f e s s i o n a l

WMROD WALL MOUNTED HEIGHT ROD KIT



**ISO13485:2003
CERTIFIED**

User Instructions

Thank you for your purchase of this Health o meter® Professional product. Please read this manual carefully, and keep it for easy reference or training.

TABLE OF CONTENTS

Caution and Warning.....	3
Specifications	3
Certifications / Disposal / Service	4
WMROD Wall Mount Instructions.....	5
WMROD Height Measurement Instructions	9
Maintenance.....	10
Warranty.....	11

Please register your product for warranty coverage at:

<http://www.homscales.com/customer/registration.aspx>

For User Instructions updates and revisions please go to:

<http://www.homscales.com>

CAUTION AND WARNING

To prevent injury and damage to your product, please follow these instructions very carefully.

- Do not transport the WMROD while headpiece or height rod is extended.
- Ensure the WMROD is secured on a flat, stable wall surface.
- Ensure that the WMROD is secured for cleaning or for reassembly.
- For accurate measurement, the WMROD must be placed on a stable wall surface.
- For accurate measurement, verify proper operation according to the procedure described in the manual before each use.
- Do not use in the presence of flammable materials.
- When not in use, ensure headpiece is at its lowest position, and is folded down against the height rod.
- When wall mounting the WMROD, ensure that the mounting location is free of electrical lines, gas pipes, and other safety hazards. Ensure that wall is of construction that can support the WMROD.
- Do not mount WMROD on a door.

SPECIFICATIONS FOR THIS PRODUCT

General


The WMROD Wall Mounted Height Rod Kit has a metal height rod with improved strength and durability. Each WMROD features a telescoping – style design that provides accurate, reliable, and repeatable height measurements, and includes features that make the process of measuring a patient's height simple, fast, and convenient. Each WMROD is set up to determine the height of a stationary patient. Height is displayed in inches (in) and centimeters (cm).

WMROD Product Specifications

Height Range	WMROD: Max 84" (213cm) / Min 23-5/8" (60cm)
Environmental	Operating temperatures: 5°C to 40°C (41°F to 104°F) Storage temperatures: 0°C to 50°C (32°F to 122°F) Atmospheric Pressure: 500hPa to 1060hPa Maximum Humidity: 85% RH in use or in storage
Physical Dimensions	WMROD: Width: 1-5/8" (4cm) Depth: 3/4" (2cm) Height: 80-3/4" (205cm)

CERTIFICATIONS / DISPOSAL / SERVICE

Certification Descriptions

	<p>Pelstar, LLC has been officially certified as the manufacturer of Health o meter® Professional medical devices. The Pelstar, LLC quality assurance system covers the development, production, sales and service of medical scales and measuring systems.</p>
---	---

WMROD Disposal



This Health o meter® Professional product must be disposed of properly. Follow the national, regional or local regulations which apply to you for disposal of domestic waste.

SERVICE INFORMATION

The WMROD Wall Mounted Height Rod Kit is shipped assembled in one carton. Carefully inspect the carton for shipping damage before unpacking. If damage is found, contact Customer Service immediately.

For further information or telephone support, please contact a Health o meter® Professional Customer Service representative at 1-800-815-6615.

WMROD WALL MOUNT INSTRUCTIONS

To avoid damaging the WMROD when unpacking it, **do not use a box cutter, knife, scissors or any sharp object when opening the protective inner packaging.** If damage is found, contact your shipper or a Health o meter® Professional representative immediately at 1-800-815-6615. Claims must be filed with the shipper as soon as possible after receipt of the package. The following information details what you will find inside the carton as you unpack the parts for assembly. Set the carton aside, for storage.

Parts List

Carton

- (2) Wall mounting brackets
- (2) Drywall anchors – plastic
- (1) Wrench
- (2) # 4 Round headed screws
- (2) Hex screws
- (1) Height rod
- (1) User instructions

Tools Required (not included)

Drill, 5/64”(1.90mm) or 17/64”(6.50mm) drill bit, Tape measure, Phillips head screwdriver, Flat head screwdriver, Hammer, 4' (125cm) bubble level.




Figure 1

WMROD WALL MOUNT INSTRUCTIONS (CONTINUED)

Wall Mounting Instructions (Refer to Figure 2 - 3)

Caution: When wall mounting the WMROD, ensure that the mounting location is free of electrical lines, gas pipes, and other safety hazards. Ensure that wall is of construction that can support the WMROD.

Note: Do not mount WMROD on a door.

Mounting on Plasterboard or Drywall (Refer to Figures 2 - 3)

Note: For hollow, drywall construction use 17/64"(6.50mm) drill bit and plastic anchors.

1. Carefully remove the product from the carton and place on a flat and dry surface.
2. Determine the location of where the height rod will be mounted.

Note: Ensure mounting area is clear of baseboard, carpeting materials, and any other obstructions.

3. Use the 17/64"(6.50mm) drill bit to drill two holes. The first hole must be positioned 6 3/16" (15.5cm) up from the floor.
4. The second hole must be centered directly above the first hole, and drilled 44 7/8" (114cm) up from the floor (see Figure 2). Use a level to position the screw hole.
5. Insert the plastic anchors into the hole. Using the two #4 round head screws, attach the wall mounting brackets to the wall (see Figures 2 & 3).
6. Screw the Hex screws into the top and bottom wall mounting brackets. Leave the Hex screw head extended at least 1/4" (7mm) for mounting the height rod.
7. Slip the height rod's key-hole mounts over the Hex screws, and ensure that the mount is secured over the Hex screw (see Figure 3).
8. Use the small end of the wrench to tighten the screws until the height rod is held firmly in place.
9. Installation is now complete and ready for use.

WMROD WALL MOUNT INSTRUCTIONS (CONTINUED)

Wall Mounting Instructions (Refer to Figure 2 - 3)

Caution: When wall mounting the WMROD, ensure that the mounting location is free of electrical lines, gas pipes, and other safety hazards. Ensure that wall is of construction that can support the WMROD.

Note: Do not mount WMROD on a door.

Mounting on Wood, Wood studs, or Metal (Refer to Figures 2 - 3)

Note: For wood, wood studs, or metal walls use 5/64" bit only

1. Carefully remove the product from the carton and place on a flat and dry surface.
2. Determine the location of where the height rod will be mounted.

Note: Ensure mounting area is clear of baseboard, carpeting materials, and any other obstructions.

3. Use the 5/64"(1.90mm) drill bit to drill two holes. The first hole must be positioned 6 3/16" (15.5cm) up from the floor.
4. The second hole must be centered directly above the first hole, and drilled 44 7/8" (114cm) up from the floor (see Figure 2). Use a level to position the screw hole.
5. Insert the plastic anchors if necessary. Using the two #4 round head screws, attach the wall mounting brackets to the wall (see Figures 2 & 3).
6. Screw the Hex screws into the top and bottom wall mounting brackets. Leave the Hex screw head extended at least 1/4" (0.7cm) for mounting the height rod.
7. Slip the height rod's key-hole mounts over the Hex screws, and ensure that the mount is secured over the Hex screw (see Figure 3).
8. Use the small end of the wrench to tighten the screws until the height rod is held firmly in place.
9. Installation is now complete and ready for use.

WMROD WALL MOUNT INSTRUCTIONS (CONTINUED)




Figure 2




Figure 3

HEIGHT MEASUREMENT INSTRUCTIONS

Height Measurements (Refer to figures 4-7)

Note: A patient needs to be a minimum of 600mm (23-5/8") in height for the height to be accurately measured.

1. Raise the headpiece above the patient's head.
2. Position the patient against the WMROD.
3. Ensure that the patient's posture is upright, and that the patient's head is straight and level.
4. Push the headpiece down until it rests on the patient's head, and tighten the black knob.
5. Determine the patient's height by reading the value on the measuring strip nearest to the line on the indicator. The indicator is located on the center top of the height rod.



Figure 4



Figure 5



Figure 6



Figure 7

MAINTENANCE

Maintenance

The following pages provide instructions for maintenance and cleaning for the WMROD product. Maintenance operations other than those described in this manual should be performed by qualified service personnel.

Caution: Before first use, or after long periods of non-use, check the product for proper operation and function. If the product does not operate correctly, refer to qualified service personnel.

1. Check overall appearance of the total product for any obvious damage, wear, and tear.

Cleaning

Proper care and cleaning is essential to ensure a long life of accurate and effective operation.

1. Clean all external surfaces with a clean damp cloth or tissue. Mild antimicrobial or antibacterial soap and water solution may be used. Dry with a clean soft cloth.
2. Do not use Isopropyl Alcohol or other solutions to clean the height rod.
3. Do not use abrasive cleaners.

WARRANTY

Limited Warranty

What does the Warranty Cover?

Health o meter® Professional product are warranted from date of purchase against defects of materials or in workmanship for a period of one (1) years. If product fails to function properly, return the product, freight prepaid and properly packed to Pelstar, LLC (see "To Get Warranty Service", below, for instructions). If the manufacturer determines that a defect of material or in workmanship exists, the customer's sole remedy will be repair or replacement of the product at no charge. Replacement will be made with a new or remanufactured product or component. If the product is no longer available, replacement may be made with a similar product of equal or greater value. All parts including repaired and replaced parts are covered only for the original warranty period.

Who is Covered?

The original purchaser of the product must have proof of purchase to receive warranty service. Please save your invoice or receipt. Pelstar dealers or retail stores selling Pelstar products do not have the right to alter, or modify or any way change the terms and conditions of this warranty.

What is Excluded?

Your warranty does not cover normal wear of parts or damage resulting from any of the following: negligent use or misuse of the product, use on improper voltage or current, use contrary to the operating instructions, abuse including tampering, damage in transit, or unauthorized repair or alternations. Further, the warranty does not cover natural disasters, such as fire, flood, hurricanes and tornadoes. This warranty gives you specific legal rights, and you may also have other rights that vary from country to country, state to state, province to province or jurisdiction to jurisdiction.

To get Warranty Service make sure you keep your sales receipt or document showing proof of purchase. Call (+1) 800-638-3722 or (+1) 708-377-0600 to receive a return authorization (RA) number, which must be included on the return label. Attach your proof of purchase to your defective product along with your name, address, daytime telephone number and description of the problem. Carefully package the product and send with shipping and insurance prepaid to:

Pelstar, LLC
Attention R/A# _____
Return Department
11800 South Austin Avenue
Alsip, IL 60803

If your product is not covered by warranty, or has been damaged, an estimate of repair costs or replacement costs will be provided to you for approval, prior to servicing or replacing.



PELSTAR, LLC

11800 South Austin Avenue • Alsip, IL 60803 • USA

1-800-638-3722 or 1-708-377-0600

PLEASE REGISTER YOUR PRODUCT FOR WARRANTY COVERAGE AT:
<http://www.homscales.com/customers/registration.aspx>

Health o meter® is a registered trademark of Sunbeam Products, Inc. used under license.
Health o meter® Professional products are manufactured, designed, and owned by Pelstar, LLC.
We reserve the right to improve, enhance, or modify Health o meter® Professional product features or specifications without notice.

© Pelstar, LLC 2012

