

HOW TO USE PRE-TARE

For Health o meter® Professional Wheelchair Scales

Know the wheelchair weight? Enter it before weighing the patient.

The Pre-Tare function allows you to enter a known weight so the scale automatically deducts it and displays only the patient's weight.

1 TURN ON THE SCALE

Make sure there is nothing on the weighing platform. Press  to turn the scale on.

2 ENTER THE PRE-TARE WEIGHT

Press and hold **TARE** until “Pre-Tare” appears at the top of the display. Use the arrow keys to enter the known weight of the wheelchair.

3 SAVE THE WEIGHT

Press **ENTER** to save the tare weight. A negative weight will be shown on the display.

4 WEIGH THE PATIENT

Have the patient move onto the scale in their wheelchair. The scale will automatically deduct the entered tare value from the weight of the patient and display only the patient's weight.

5 RECORD AND CLEAR THE WEIGHT

Record the patient's weight. Then the patient may move off of the scale. To remove the tare and ready the scale for the next patient, press **TARE**.



Pre-Tare can be performed within the following weight limits: 0.2 lb - 500 lb / 0.1 kg – 226 kg)