



US Patent #D667323

### Advantages

- ✓ **EMR Connectivity**  
Scale has the ability to interface with a computer, EMR software or other electronic device via various Health o meter<sup>®</sup> Professional Connectivity Solutions
- ✓ **350° Swivel Display**  
Eye-level rotating display can be read from almost any angle and allows for ease of use and patient privacy
- ✓ **Digital Height Rod**  
Fast and easy-to-use digital height rod provides option to capture height measurement electronically and accurately
- ✓ **Value-Added Features**  
Operating functions (LB/KG Conversion, LB/KG Lock Out, BMI, Zero Out/Tare, Hold/Release, Auto Zero, Auto Off) included at no additional cost
- ✓ **Calculates Body Mass Index (BMI)**  
Scale functions include the ability to calculate the patient's Body Mass Index (BMI)
- ✓ **Low Profile Platform**  
Low platform height requires only minimal step up, offering patient a more comfortable and stable weighing experience

# Heavy Duty Eye Level Digital Scale with Integral Digital Height Rod

600KL

**Health o meter®**  
P r o f e s s i o n a l



Model	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
600KL	52" x 26.375" x 6.375" / 5.06	38.25	854896003115



US Patent #D667323

## Specifications

- Capacity: 600 lb / 272 kg
- Resolution: 0.2 lb / 0.1 kg
- EMR connectivity via USB
- Includes Connectivity Kit
- 1" LCD display, 35° left / right swivel (350° with height rod fully retracted) for ease of viewing
- Functions: LB / KG Conversion, LB /KG Lock Out, Body Mass Index (BMI), Zero Out / Tare, Hold / Release, Auto Zero, Auto Off
- Platform Size: 14 1/4" (w) x 14 1/4" (d) x 2 5/8" (h) / 362 mm (w) x 362 mm (d) x 67 mm (h)
- Digital Height Rod: 24" - 88" / 61 cm - 223 cm
- (2) Wheels for ease of movement
- 120V adapter included, 6 C-cell batteries not included
- 2 Year Limited Warranty