

Powered by

**VALHALLA**  
SCIENTIFIC

Body Composition Scale

**ADVANTAGES**

✔ **15 Second Measurement Time**

On screen display delivers rapid access to weight, body fat %, body water %, lean mass, Basal Metabolic Rate, & Body Mass Index. A full illustrated printout can be generated in a minute

✔ **Large Platform**

The spacious platform allows patients to find a comfortable and secure stance.

✔ **Few Pretest Requirements**

Analysis can be performed at any time of day and under almost any circumstances, resulting in the greatest flexibility of any body composition analyzer.

✔ **High Capacity**

The scale features an industry leading 1000 lb capacity, that provides the flexibility and accuracy to weigh any patient.

✔ **Patient Report Templates\***

The G6 Series scales offer users 4 graphic printout formats to use with their patients: Adult/Bariatric, Adult Spanish, Pediatric and Fitness.

\* User must connect scale to a printer to provide personalized graphic printout. Printer is optional and sold separately, not available through Health o meter® Professional scales.

 **ScaleSurance Certified**

Extended Warranty Program that offers an extra 2 years of coverage

Powered by

**VALHALLA**  
SCIENTIFIC

Body Composition Scale

Model	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
BCS-G6	46" x 19" x 10" / 10	37 lb	854896003702

**Specifications**

- Capacity: 1000 lb / 453.5 Kg
- Resolution: 0.1 lb / 0.1 kg
- Platform Size: 14" (w) x 16" (d) x 1.875" (h)
- Display: 6-window display (1/2" digits) includes weight, % body fat, % body water, % muscle mass, BMI, REE
- Power Source: 110-220V power adapter (included)
- Optional printer, sold separately, not available through Health o meter® Professional scales
- 1 Year Limited Warranty - Eligible for ScaleSurance Extended Warranty

**Bio-Impedance Specifications**

- Electrode Method: Tetrapolar System
- Frequency: 50kHz
- Test Current: 500µA
- Conversion Rate: 750msec



<b>BCS-G6-adult</b> Body Comp Scale	The whole-body, adult only scale is ideal for practices that exclusively see patients 18 and over	854896003702
<b>BCS-G6-limbs</b> Body Comp Scale	The adult segmented version provides patients with body composition information about their torso as well as specific information about the composition of their arms and legs	854896003719
<b>BCS-G6-DUO</b> Body Comp Scale	The whole-body adult/pediatric version is a good choice for providers who see adults and children	854896003733
<b>BCS-G61-upper</b> Body Comp Scale	The adult upper-body-only configuration provides a partial body composition analysis for patients who do not want to remove their shoes and socks	854896003726
<b>IPO-ADU</b>	Patient Report Template - Adult (5-page report, 25 reports/package)	854896003740
<b>IPO-PED</b>	Patient Report Template - Pediatric (7-page report, 25 reports/package)	854896003757
<b>IPO-ATH</b>	Patient Report Template - Athletic (3-page report, 50 reports/package)	854896003764
<b>IPO-ADU-SP</b>	Spanish Language Patient Report Template - Adult (5-page report, 25 reports/package)	854896003948
<b>BCS-BROCHURE</b>	Waiting Room Patient Brochure	