

THE 349KLX

DIGITAL FLOOR SCALE WITH REMOTE DISPLAY

The 349KLX is a remote display digital scale that supports up to 400 pounds and provides the optimum in quality, reliability, and ease of use. The remote display provides flexibility with the ability to be mounted on a wall or placed on a table. It can be easily moved from location to location, to adapt to various workplace environments.

Designed to meet the needs of doctor's offices, hospitals, and clinics, the scale's large easy-to-read display, and the large non-skid low profile platform requires a minimal step-up, providing a safe weighing environment. For easy portability and storage, there is an optional carrying case.

The user-friendly digital interface includes features such as Zero, Tare, Hold/ Release, Auto-off and more. The exclusive Everlock® function allows the measuring unit (LB or KG) to be permanently locked on the scale.



THERE'S STRENGTH IN OUR NUMBERS.

When you add up the numbers, it's not surprising that the Health o meter® Professional Scales is the industry-leading name in medical scales.

120+ LOCAL
FIELD SPECIALISTS

9' REMOTE
CABLE LENGTH

#1 TRUSTED
BRAND

100 YEARS
OF SERVICE

Health o meter®

Professional

Health o meter® Professional 349KLX

weigh**easier**

TECHNICAL SPECIFICATIONS

Capacity

400 lb / 180 kg

Resolution

0.2 lb / 0.1 kg

Display

1" / 25 mm, LCD display

Power Source

6 AA batteries included,
optional adapter
(order # ADPT40)

Functions

LB/KG conversion
LB/KG lock out
Zero
Tare
Hold/release
Reweigh
Auto zero
Auto off
Variable auto off time
Parent/child weighing function
Enabling unit view with unit of measure locked
Enable/disable backlight display with 4 levels of brightness
Everlock®

Platform Size (w x d x h)

12 1/2" x 12" x 1 7/8"
(318 mm x 305 mm x 48 mm)

Product Footprint - Base (w x d x h)

12 1/2" x 12" " x 1 7/8"
(318 mm x 305 mm x 48 mm)

Product Weight

10 lb / 4.5 kg

Portability

64771 - optional carrying case

Remote Cable Length

9 ft / 2743 mm

Warranty

2-year limited manufacturer warranty

ORDERING INFORMATION

MODEL	DESCRIPTION
349KLX	Digital floor scale with remote display, LB & KG
349KGX	Digital floor scale with remote display, KG only
SS-349KLX	ScaleSurance® 2 year extended warranty for 349KLX series*
349KLXAD	349KLX with adapter ADPT40 included
64771	Carrying case for 349KLX series of scales
ADPT40	Optional power adapter for 349KLX series of scales
PLUGSET30	Plug adapters for use outside North America



64771

Carrying Case for 349KLX Series of Scales



To learn more, call **1-800-815-6615**
or visit **hom scales.com/349KLX**

Health o meter®
Professional

*Not available in all countries.