

THE PTS-1000KG-C

PATIENT TRANSFER SCALE, KG ONLY, CANADA ONLY

The PTS-1000KG-C is a transfer board with a built-in medical scale that provides a weight measurement when transferring a patient. A kilograms only Patient Transfer Scale offers a safe solution for facilities standardizing on metric only measurements. All departments can benefit from the Patient Transfer Scale, particularly Stroke Units, the Emergency Room, ICU, and Radiology.

Designed by Clinicians for Clinicians the Patient Transfer Scale weighs critical or immobile patients instantly and accurately helping shorten door-to-needle times. A kilograms only Patient Transfer Scale is especially useful for Stroke Units where quick and accurate patient weights are crucial and where dosages of critical medications are based on kilogram measurements. It offers a superior alternative to other weighing equipment, such as in-bed scales, that are often not calibrated or inaccurate.

The scale is light and portable and can be utilized as part of an existing patient transfer process when a weight is needed, and will not disrupt the patient flow process.



THERE'S STRENGTH IN OUR NUMBERS.

When you add up the numbers, it's not surprising that Health o meter® Professional Scales is the industry-leading name in medical scales.

120+ LOCAL
FIELD SPECIALISTS

#1 TRUSTED
BRAND

100 YEARS
OF SERVICE



Health o meter®
Professional

Health o meter® Professional PTS-1000KG-C

weigh**easier**®

TECHNICAL SPECIFICATIONS

Capacity

250 kg

Resolution

0.5 kg

Display

3" x 1" / 76 mm x 25 mm Backlight LCD

Power Source

3.7V rechargeable battery

Functions

Zero
Hold/release
Auto off

Weighing Mat Size (w x d)

71" x 21"
(1805mm x 533mm)

Product Footprint (w x d x h)

71" x 27 1/2" x 1"
(1805mm x 700mm x 30mm)

Product Weight

25 lb / 11.4 kg

Measuring Tape

Height guide shown on mat 100cm – 180cm
shown at 5 cm

Warranty

2-year limited manufacturer warranty

ORDERING INFORMATION

MODEL	DESCRIPTION
PTS-1000KG-C*	Patient Transfer Scale, KG only, Canada only
SS-PTS-1000KL**	ScaleSurance® 2 year extended warranty for Patient Transfer Scale



Always use a Slide Sheet
with the Patient Transfer Scale!



Check out the Single Patient Use
and Reusable Slide Sheets.

To learn more, call **1-800-815-6615**
or visit [homcales.com/PTS-1000KG-C](https://www.homcales.com/PTS-1000KG-C)

Health o meter®
Professional

* Evaluated to SPE 3000 for sale in Canada
** Not available in all countries.