



The 499KGHR is a compact version of our most popular scale, the 500KG. The waist-high display is easy to read and the scale is an attractive, robust option for use in medical offices. With a height rod included, the 499KGHR is an economical choice for almost any environment.

Advantages

- ✓ **EMR Connectivity**
Scale has the ability to interface with a computer, EMR software or other electronic device via various Health o meter[®] Professional Connectivity Solutions
- ✓ **Calculates Body Mass Index (BMI)**
Scale functions include the ability to calculate the patient's BMI
- ✓ **Low Profile Platform**
Low platform height requires only minimal step up, offering patient a more comfortable and stable weighing experience
- ✓ **Value-Added Features**
Operating functions (BMI, Zero, Tare, Hold/Release, Auto Zero, Auto Off) included at no additional cost
- ✓ **Compact Size**
The slim design of the scale and waist-high display require a limited amount of space



Specifications

- Capacity: 220 kg
- Resolution: 0.1 kg
- EMR connectivity via USB, optional Pelstar® wireless technology
- 1" / 25 mm, LCD display
- Functions: BMI, Zero, Tare, Hold / Release, Auto Zero, Auto Off
- Platform Size (w x d x h):
13 ¾" x 16 ½" x 2 ¾" / 349 mm x 419 mm x 60 mm
- Product weight: 21 lb / 9 kg
- Height rod: 11 7/8" - 84" / 30 cm - 213 cm
- Height rod graduation: 1/8" / 1 mm
- 2 wheels for ease of movement
- 6 AA batteries included, optional 120V adapter (order # ADPT31)
- 2 year limited warranty

Model	Description
499KGHR	Waist high digital scale, kilograms only, height rod included
499KLHR	Waist high digital scale, pounds & kilograms, height rod included
499KGHRAD	499KG with height rod & adapter ADPT31 included
499KLHRAD	499KL with height rod & adapter ADPT31 included
SS-499KL	ScaleSurance 2 year extended warranty for 499 scale models*
ADPT31	Optional 120V adapter
PLUGSET30	Plug adapters for use outside North America
PELSTARONE-P	Wireless thermal printer