



US Patent #D667323

Improve workflow with this top-of-the-line physician office scale. The built-in digital height rod eliminates the need to manually enter height for computing BMI. The 350° swivel head allows for line-of-sight viewing and patient privacy.

Advantages

- ✓ **EMR Capable**
Scale interfaces with a computer, EMR software or other electronic device via Health o meter® Professional Connectivity Solutions
- ✓ **350° Swivel Display**
Eye-level rotating display can be read from almost any angle and allows for ease of use and patient privacy
- ✓ **Built-in Digital Height Rod**
The easy-to-use, state of the art digital height rod provides measurements from 23 5/8" - 88" / 60 cm - 222 cm
- ✓ **Value-Added Features**
Operating functions (BMI, Zero, Tare, Hold/Release, Auto Zero, Auto Off) included at no additional cost
- ✓ **Calculates Body Mass Index (BMI)**
Scale functions include the ability to calculate the patient's BMI
- ✓ **Patented Design**
The 600KG's patented design provides a user-friendly and intuitive interface for measuring patient weight, height and BMI with a convenient rotating display head



US Patent #D667323

Specifications

- Capacity: 272 kg
- Resolution: 0.1 kg
- Platform Size (w x d x h):
14 ¼" x 14 ¼" x 2 ⅝"
362 mm x 362 mm x 67 mm
- Product Footprint (w x h x d):
14 ¾" x 24 ⅛" x 54 ⅜"
375 mm x 613 mm x 1381 mm
- Product Weight: 28 lb / 13 kg
- Digital Height Rod: 23 ⅝" - 88" / 60 cm - 222 cm
- Height Rod Resolution: 0.1" / 0.1cm
- Connectivity: USB, optional Pelstar wireless technology
- Display: 1" / 25 mm LCD, 35° left / right swivel
(350° with height rod fully retracted) for ease of viewing
- Functions: BMI, Zero, Tare, Hold / Release, Auto Zero, Auto Off
- 2 wheels for ease of movement
- 100 - 240V adapter (ADPT50) included,
6 C-cell batteries (not included)
- 2 year limited warranty

Model	Description
600KG-BT	600KG with built-in Pelstar wireless technology
SS-600KL	ScaleSurance 2 year extended warranty for 600KG*
SS-600KL-BT	ScaleSurance 2 year extended warranty for 600KG-BT*
ADPT50	Adapter included, item # for replacement or spare
PLUGSET50	Plug adapters for use outside North America
PELSTARONE-P	Wireless thermal printer