

# Health o meter®


P r o f e s s i o n a l

## Wall Mounted Height Rod

PORTROD



- Provides accurate measurement of a patient's height
- Large headpiece with clear indicator
- Height is displayed in centimeters (cm), inches (in) and feet/inches
- Range: Max. 213cm (83") / Min. 60cm (24")
- Accuracy: +/- 3mm (1/8")
- Physical Dimensions: Width: 7cm (2.75")
- Depth: 1.3cm (.5")
- Height: 152.4cm (60")



Health o meter® Professional Scales  
9500 West 55th Street USA  
McCook, IL 60525-7110  
Phone: 1-800-815-6615  
Fax: 1-800-654-7330

**Health o meter®**  
Professional  
[www.homscales.com](http://www.homscales.com)