

THE 498KL

REMOTE DISPLAY DIGITAL FLOOR SCALE

The 498KL is a remote display digital scale that supports up to 660 pounds and provides the optimum in quality, reliability, and ease of use. The remote display provides flexibility with the ability to be mounted on a wall or placed on a table. It can be easily moved from location to location, to adapt to various workplace environments.

Designed to meet the needs of doctor's offices, hospitals, and clinics, the scale's large easy-to-read display, and the large non-skid low profile platform requires a minimal step-up, providing a safe weighing environment. For easy portability and storage, there is an optional carrying case.

The user-friendly digital interface includes features such as Zero, Tare, Hold/Release, Auto-off and more. The exclusive Everlock® function allows the measuring unit (LB or KG) to be permanently locked on the scale.

The 498KL provides easy integration and seamless communication with today's EMR applications with the option to incorporate wireless connectivity for interfacing with computers, EMR software, or other electronic devices.



THERE'S STRENGTH IN OUR NUMBERS.

When you add up the numbers, it's not surprising that the Health o meter® Professional Scales is the industry-leading name in medical scales.

120+ LOCAL
FIELD SPECIALISTS

9' REMOTE
CABLE LENGTH

660 LB
CAPACITY

#1 TRUSTED
BRAND

100 YEARS
OF SERVICE

Health o meter®
Professional

Health o meter® Professional 498KL

weigh**easier**

TECHNICAL SPECIFICATIONS

Capacity

660 lb / 300 kg

Resolution

0.2 lb / 0.1 kg

EMR Connectivity

Via optional Pelstar® wireless technology

Display

1" / 25 mm, LCD display

Power Source

6 AA batteries included, optional adapter (order # ADPT31)

Functions

LB/KG conversion
LB/KG lock out
BMI
Zero
Tare
Pre-tare
Hold/release
Reweigh
Auto zero
Auto off
Variable auto off time
Audible/mute sound option
Parent/child weighing function
Enabling unit view with unit of measure locked
Everlock®

Platform Size (w x d x h)

13 3/4" x 16 1/2" x 2 3/8"
(349 mm x 419 mm x 60 mm)

Product Footprint - Base (w x d x h)

13 3/4" x 16 1/2" x 2 3/8"
(349 mm x 419 mm x 60 mm)

Product Weight

14 lb / 6.4 kg

Portability

64771 - optional carrying case

Remote Cable Length

9 ft / 2743 mm

Warranty

2-year limited manufacturer warranty

ORDERING INFORMATION

MODEL	DESCRIPTION
498KL	Remote display digital floor scale, LB & KG
498KG	Remote display digital floor scale, KG only
498KLAD	498KL with adapter ADPT31 included
SS-498KL	ScaleSurance® 2 year extended warranty for 498KL series*
64771	Carrying case for 498KL series of scales
ADPT31	Optional power adapter for 498KL series of scales
PLUGSET30	Plug adapters for use outside North America
BTPRINTERKIT-01	Wireless thermal printer kit

Scale also includes connectivity capabilities. To learn about connectivity options, visit homscal.es.com/ConnectivitySolutions



64771
Carrying Case for 498KL Series of Scales

To learn more, call **1-800-815-6615**
or visit homscal.es.com/498KL



*Not available in all countries.

Health o meter®
Professional