

Remote Display can be mounted on a wall or placed on a table

Advantages

Calculates Body Mass Index (BMI)
Scale functions include the ability to
calculate the patient's BMI

Low Profile Platform

Low platform height and wide weighing surface offer patient a more comfortable and stable weighing experience **✓ Value-Added Features**

Operating functions (LB/KG Conversion, LB/KG Lock Out, BMI, Zero Out/Tare, Hold/Release, Auto Zero, Auto Off, Parent/Child Weighing Function) included at no additional cost

Remote Display

The remote display provides versatility with the ability to be mounted on a wall or placed on a table





Remote Display can be mounted on a wall or placed on a table

Model	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
498KL	25.875" x 17.375" x 5.5" / 1.43	16.25 lb	892076002718



Specifications

- Capacity: 500 lb / 220 kg
- Resolution: 0.2 lb / 0.1 kg
- Platform Size: 13 ¾" (w) x 16 ½" (d) x 2 ¾" (h) / 349 mm (w) x 419 mm (d) x 60 mm (h)
- Display: 1" LCD
- Optional Carrying Case: Item# 64771
- Power Source: 6 AA batteries included, Optional adapter (order # ADPT31)
- Functions: LB / KG Conversion, LB / KG Lock Out, Body Mass Index (BMI), Zero Out / Tare, Hold / Release, Auto Zero, Auto Off, Parent / Child Weighing Function
- 2 Year Limited Warranty



Carrying Case for 498KL Item# 64771



Health o meter® Professional Scales 9500 West 55th Street USA McCook, IL 60525-7110 Phone: 1-800-815-6615 Fax: 1-800-654-7330