



Advantages

- ✓ **125 kg Capacity Floor Scale**
The high capacity of these floor scales provide more versatile use
- ✓ **Easy-to-read Display**
Large, easy-to-read display is conveniently located at the top of the scale
- ✓ **Small-Step Platform**
Low platform height requires only minimal step up, offering patient a more comfortable and stable weighing experience



Specifications

- Capacity: 125 kg
- Graduation: 0.5 kg
- Dial Size: 2" / 51 mm
- Platform Size (w x d x h):
10 1/4" x 10" x 1 7/8" /
260 mm x 254 mm x 48 mm
- Product Footprint (w x d x h):
10 1/4" x 10" x 2" /
260 mm x 254 mm x 51 mm
- Product Weight: 2 lb / 0.9 kg
- Optional Carrying Case Available (order # 64771)
- 1 Year Limited Warranty



Carrying Case for 100KG
Item# 64771