THE 1100KL-EHR

DIGITAL PLATFORM SCALE WITH EXTRA WIDE HANDRAILS AND DIGITAL HEIGHT ROD

The all-digital 1100KL-EHR combines the industry-leading bariatric scale that supports up to 1000 pounds and provides the optimum quality, reliability, and ease of use with a state-of-the-art digital height rod. The market-leading height rod provides the most accurate height measurement in its class, with a graduation/resolution of 1 mm/0.1" and a full measuring range of 11 $\frac{3}{4}$ " – 79 $\frac{1}{2}$ " / 30cm – 202cm. The extra-wide "live" handrails allow the patient to hold on for stability while the scale accurately calculates weight and height.

Built-In Reset Feature

Proprietary technology eliminates mechanical sensors, reducing failures, increasing accuracy and reliability

Designed to meet the needs of doctor's offices, hospitals, clinics, and acute and longterm care facilities, the scale's keypad features antimicrobial protection, a display head that swivels 180°, a large non-skid low profile platform, and wheels that enable portability

The scale's high-contrast color screen provides a user-friendly interface that includes BMI, hold/release, reweigh, recall, zero, KG/LB, tare, pre-tare, and more. The exclusive Everlock® function allows the measuring unit (LB or KG) to be permanently locked on the scale.

The scale provides easy integration and seamless communication with today's EMR applications. The 1100KL-EHR includes a USB interface for connectivity with a computer, EMR software, or other electronic devices via various Health o meter[®] Professional Connectivity Solutions.



Health o meter® Professional 1100KL-EHR

Industry leading bariatric scale as reported by Global Industry Analysts, Inc. Patient Scales market survey and Global Healthcare Exchange LLC, (GHX) sales data.

THERE'S STRENGTH IN OUR NUMBERS.

When you add up the numbers on our scale, it's not surprising that the 1100KL-EHR is the industry leading bariatric scale.

120+ LOCAL FIELD SPECIALISTS

100% Connectivity

#1 TRUSTED BRAND

100 YEARS OF SERVICE

weigh**easier**®

TECHNICAL SPECIFICATIONS

Capacity 1000 lb / 454 kg

Resolution 0.2 lb / 0.1 kg

EMR Connectivity Via USB, optional Pelstar® wireless technology

Display 1 ½" / 38 mm High-Contrast Color TFT-LCD Screen

Power Source 100-240V adapter (ADPT30) included, 6 D batteries not included

Digital Height Range 3ft 11.2" – 6ft 7.5" / 120 cm - 202 cm Resolution: 1 mm / 0.1"

Mechanical Height Range 11 ¾" – 79 ½" / 30cm – 202cm Graduation: 1mm / 1/16 inch

Functions

LB/KG conversion LB/KG lock out BMI Zero Tare Pre-tare Hold/release Reweigh Recall Auto zero Auto off Time/date Variable auto off time Audible/mute sound option Parent/child weighing function Everlock[®] - permanent weight lock

Platform Size (w x d x h) 22" x 15 ³4" x 2 ¹⁄2" (559 mm x 400 mm x 64 mm)

Product Footprint (height rod extended) 31 ³/₄" x 33 ¹/₂" x 82 ¹/₂" (806 mm x 851 mm x 2096 mm)

Product Footprint

(height rod in home position) 31 ¾" x 33 ½" x 50 ¼" (806 mm x 851 mm x 1276 mm)

Product Weight 58 lb / 26 kg

Portability 2 wheels for ease of movement

Warranty 2-year limited manufacturer warranty

The most accurate digital height rod in its class with a graduation/ resolution of 1mm / 0.1"

Industry-Leading Accuracy

ORDERING INFORMATION

MODEL	DESCRIPTION
1100KL-EHR	Digital platform scale with extra wide handrails and digital height rod, LB & KG
1100KL-BT-EHR	1100KL-EHR with Pelstar® wireless technology
SS-1100KL-EHR	ScaleSurance® 2 year extended warranty
ADPT30	Adapter included, item # for replacement or spare
PLUGSET30	Plug adapters for use outside North America
PELSTARONE-P	Wireless thermal printer

Height Rod Level to Platform

Bottom of height rod at the platform level, ensuring caregiver confidence and accurate measurement every time





