

THE 1100KL-EHR

DIGITAL PLATFORM SCALE WITH EXTRA WIDE HANDRAILS AND DIGITAL HEIGHT ROD

The Certified ADA Compliant 1100KL-EHR combines the industry-leading bariatric scale that supports up to 1000 pounds and provides the optimum quality, reliability, and ease of use with a state-of-the-art digital height rod. The market-leading height rod provides the most accurate height measurement in its class, with a graduation/resolution of 1mm/0.1" and a full measuring range of 11 ¾" – 79 ½" / 30cm – 202cm. The extra-wide "live" handrails on both sides of the platform allow the patient to hold on for stability while the scale accurately calculates weight and height.



Built-In Reset Feature

Proprietary technology eliminates mechanical sensors, reducing failures, increasing accuracy and reliability

Designed to meet the needs of doctor's offices, hospitals, clinics, and acute and long-term care facilities, the scale's keypad features antimicrobial protection, a display head that swivels 180°, a large slip-resistant ow-profile platform, and wheels that enable portability.

The scale's high-contrast color screen provides a user-friendly interface that includes BMI, hold/release, reweigh, recall, zero, KG/LB, tare, pre-tare, and more. The exclusive Everlock® function allows the measuring unit (LB or KG) to be permanently locked on the scale. The scale provides easy integration and seamless communication with today's EMR applications, the 1100KL-EHR includes a USB interface for connectivity with a computer, EMR software, or other electronic devices via various Health o meter® Professional Connectivity Solutions.

The 1100KL-EHR digital platform scale with extra wide handrails and digital height rod was designed to meet the US Access Boards Standards: Medical facilities operated by state and local governments, also known as public entities, must comply with the MDE Standards for Accessible Medical Diagnostic Equipment, which is enforced by the Department of Justice.

- Slip-resistant platform
- Handrail supports on both sides of platform
- 34" – 38" max. handrail height from platform

Reviewed for ADA compliance and
Platinum certified by Corada

Health o meter®
Professional

Health o meter® Professional 1100KL-EHR

Industry leading bariatric scale as reported by Global Industry Analysts, Inc.
Patient Scales market survey and Global Healthcare Exchange LLC, (GHX) sales data.



REVIEWED FOR ADA COMPLIANCE
CERTIFICATION N° 5776764

THERE'S STRENGTH IN OUR NUMBERS.

When you add up the numbers on our scale, it's not surprising that the 1100KL-EHR is the industry leading bariatric scale.

**100% CERTIFIED
FOR ADA COMPLIANCE**

**120+ LOCAL
FIELD SPECIALISTS**

**100%
CONNECTIVITY**

**#1 TRUSTED
BRAND**

**100 YEARS
OF SERVICE**

weigh**easier**®



TECHNICAL SPECIFICATIONS

Capacity

1000 lb / 454 kg

Resolution

0.2 lb / 0.1 kg

EMR Connectivity

Via USB, optional Pelstar® wireless technology

Display

Digit Size: 3/8" / 10 mm
Screen Size: 2 3/4" x 1 1/2" High-Contrast
Color TFT-LCD Screen

Power Source

100-240V adapter (ADPT30)
included, 6 D batteries not included

Digital Height Range

3ft 11.2" – 6ft 7.5" / 120 cm - 202 cm
Resolution: 1 mm / 0.1"

Mechanical Height Range

11 3/4" – 79 1/2" / 30cm – 202cm
Graduation: 1 mm / 1/16 inch

Functions

LB/KG conversion
LB/KG lock out
BMI
Zero
Tare
Pre-tare
Hold/release
Reweigh
Recall
Auto zero
Auto off
Time/date
Variable auto off time
Audible/mute sound option
Parent/child weighing function
Everlock® - permanent weight lock

Platform Size (w x d x h)

22" x 15 3/4" x 2 1/2"
(559 mm x 400 mm x 64 mm)

Product Footprint

(height rod extended)
31 3/4" x 33 1/2" x 82 1/2"
(806 mm x 851 mm x 2096 mm)

Product Footprint

(height rod in home position)
31 3/4" x 33 1/2" x 50 1/4"
(806 mm x 851 mm x 1276 mm)

Product Weight

58 lb / 26 kg

Portability

2 wheels for ease of movement

Warranty

2-year limited manufacturer warranty

ORDERING INFORMATION

MODEL	DESCRIPTION
1100KL-EHR	Digital platform scale with extra wide handrails and digital height rod, LB & KG
1100KL-BT-EHR	1100KL-EHR with Pelstar® wireless technology
SS-1100KL-EHR	ScaleSurance® 2 year extended warranty*
ADPT30	Adapter included, item # for replacement or spare
PLUGSET30	Plug adapters for use outside North America
BTPRINTERKIT-01	Wireless thermal printer kit

Scale also includes connectivity capabilities. To learn about connectivity options, visit hom scales.com/ConnectivitySolutions

Height Rod Level to Platform

Bottom of height rod at the platform level, ensuring caregiver confidence and accurate measurement every time

Industry-Leading Accuracy

The most accurate digital height rod in its class with a graduation/resolution of 1 mm / 0.1"



To learn more, call **1-800-815-6615**
or visit hom scales.com/1100KL-EHR

Health o meter®
Professional



*Not available in all countries.

1100KL-EHR_F&B_20250825