

THE 1100KL

DIGITAL PLATFORM SCALE WITH EXTRA WIDE HANDRAILS

The all-digital 1100KL is the industry leading bariatric scale that supports up to 1000 pounds and provides the optimum in quality, reliability, and ease of use. The extra wide "live" handrails allow the patient to hold on for stability while the scale accurately calculates weight. Designed to meet the needs of doctor's offices, hospitals, clinics, and acute and long-term care facilities, the scale's keypad features antimicrobial protection, a display head that swivels 180 °, a large non-skid low profile platform, and wheels that enable portability.

The scale's high-contrast color screen provides a user friendly interface that includes BMI, hold/release, reweigh, recall, zero, KG/LB, tare, pre-tare and more. The exclusive Everlock® function allows the measuring unit (LB or KG) to be permanently locked on the scale. The scale provides easy integration and seamless communication with today's EMR applications, the 1100KL includes a USB interface for connectivity with a computer, EMR software, or other electronic devices via various Health o meter® Professional Connectivity Solutions.



THERE'S STRENGTH IN OUR NUMBERS.

When you add up the numbers on our scale, it's not surprising that the 1100KL is the industry leading bariatric scale.

120+ LOCAL
FIELD SPECIALISTS

100%
CONNECTIVITY

#1 TRUSTED
BRAND

100 YEARS
OF SERVICE

Health o meter®
Professional

Health o meter® Professional 1100KL

Industry leading bariatric scale as reported by Global Industry Analysts, Inc. Patient Scales market survey and Global Healthcare Exchange LLC, (GHX) sales data.

weigh**easier**®

TECHNICAL SPECIFICATIONS

Capacity

1000 lb / 454 kg

Resolution

0.2 lb / 0.1 kg

EMR connectivity

Via USB, optional Pelstar® wireless technology

Display

Digit Size: 3/8" / 10 mm

Screen Size: 2 3/4" x 1 1/2" High-Contrast

Color TFT-LCD Screen

Power Source

100-240V adapter (ADPT30)

included, 6 D batteries not included

Functions

LB/KG conversion

LB/KG lock out

BMI

Zero

Tare

Pre-tare

Hold/release

Reweigh

Recall

Auto zero

Auto off

Time/date

Variable auto off time

Audible/mute sound option

Parent/child weighing function

Everlock®

Platform Size (w x d x h)

22" x 15 3/4" x 2 1/2"

(559 mm x 400 mm x 64 mm)

Product Footprint (w x d x h)

31 3/4" x 33 1/2" x 49 3/4"

(806 mm x 851 mm x 1264 mm)

Product Weight

52 lb / 23 kg

Portability

2 wheels for ease of movement

Warranty

2-year limited manufacturer warranty

ORDERING INFORMATION

MODEL	DESCRIPTION
1100KL	Digital platform scale with extra wide handrails, LB & KG
1100KG	Digital platform scale with extra wide handrails, KG only
1100KL-BT	1100KL with built-in Pelstar® wireless technology
1100KG-BT	1100KG with built-in Pelstar® wireless technology
1100KL-EHR	1100KL with digital height rod
1100KL-BT-EHR	1100KL-EHR with built-in Pelstar® wireless technology
1100KLHR	1100KL with mechanical height rod
1100KGHR	1100KG with optional height rod
245EHR-1100	Optional digital height rod for 1100 series*
PROPLUSROD	Optional mechanical height rod for 1100 series
SS-1100KL	ScaleSurance® 2 year extended warranty for 1100 series**
ADPT30	Adapter included, item # for replacement or spare
PLUGSET30	Plug adapters for use outside North America
BTPRINTERKIT-01	Wireless thermal printer kit

To learn more, call **1-800-815-6615**
or visit **homscapes.com/1100KL**



* Compatible with 1100 series scales including a date code with the last two digits of 21 or higher. An upgrade option is available for models including a date code with the last two digits of 18 or higher. For assistance, contact Customer Service at 800-815-6615.

**Not available in all countries.



Health o meter®
Professional

1100KL_F&B_20231116