

## THE 1100KL

# DIGITAL PLATFORM SCALE WITH EXTRA WIDE HANDRAILS

The Certified ADA Compliant 1100KL is the industry leading bariatric scale that supports up to 1000 pounds and provides the optimum in quality, reliability, and ease of use. The extra wide "live" handrails on both sides of the platform allow the patient to hold on for stability while the scale accurately calculates weight. Designed to meet the needs of doctor's offices, hospitals, clinics, and acute and long-term care facilities, the scale's keypad features antimicrobial protection, a display head that swivels 180°, a large slip-resistant low-profile platform, and wheels that enable portability.

The scale's high-contrast color screen provides a user friendly interface that includes BMI, hold/release, reweigh, recall, zero, KG/LB, tare, pre-tare and more. The exclusive Everlock® function allows the measuring unit (LB or KG) to be permanently locked on the scale. The scale provides easy integration and seamless communication with today's EMR applications, the 1100KL includes a USB interface for connectivity with a computer, EMR software, or other electronic devices via various Health o meter® Professional Connectivity Solutions.

The 1100KL digital platform scale with extra wide handrails was designed to meet the US Access Boards Standards: Medical facilities operated by state and local governments, also known as public entities, must comply with the MDE Standards for Accessible Medical Diagnostic Equipment, which is enforced by the Department of Justice.

- Slip-resistant platform
- Handrail supports on both sides of platform
- 34" – 38" max. handrail height from platform

Reviewed for ADA compliance and  
Platinum certified by Corada

**Health o meter®**  
Professional

Health o meter® Professional 1100KL

Industry leading bariatric scale as reported by Global Industry Analysts, Inc. Patient Scales market survey and Global Healthcare Exchange LLC, (GHX) sales data.



## THERE'S STRENGTH IN OUR NUMBERS.

When you add up the numbers on our scale, it's not surprising that the 1100KL is the industry leading bariatric scale.

**100% CERTIFIED**  
FOR ADA COMPLIANCE


**120+ LOCAL**  
FIELD SPECIALISTS

**100%**  
CONNECTIVITY

**#1 TRUSTED**  
BRAND

**100 YEARS**  
OF SERVICE

weigh**easier**®



## TECHNICAL SPECIFICATIONS

### Capacity

1000 lb / 454 kg

### Resolution

0.2 lb / 0.1 kg

### EMR connectivity

Via USB, optional Pelstar® wireless technology

### Display

Digit Size: 3/8" / 10 mm  
Screen Size: 2 3/4" x 1 1/2" High-Contrast  
Color TFT-LCD Screen

### Power Source

100-240V adapter (ADPT30)  
included, 6 D batteries not included

### Functions

LB/KG conversion  
LB/KG lock out  
BMI  
Zero  
Tare  
Pre-tare  
Hold/release  
Reweigh  
Recall  
Auto zero  
Auto off  
Time/date  
Variable auto off time  
Audible/mute sound option  
Parent/child weighing function  
Everlock® - permanent unit lock

### Platform Size (w x d x h)

22" x 15 3/4" x 2 1/2"  
(559 mm x 400 mm x 64 mm)

### Product Footprint (w x d x h)

31 3/4" x 33 1/2" x 49 3/4"  
(806 mm x 851 mm x 1264 mm)

### Product Weight

52 lb / 23 kg

### Portability

2 wheels for ease of movement

### Warranty

2-year limited manufacturer warranty

## ORDERING INFORMATION

MODEL	DESCRIPTION
1100KL	Digital platform scale with extra wide handrails, LB & KG
1100KG	Digital platform scale with extra wide handrails, KG only
1100KL-BT	1100KL with built-in Pelstar® wireless technology
1100KG-BT	1100KG with built-in Pelstar® wireless technology
1100KL-EHR	1100KL with digital height rod
1100KL-BT-EHR	1100KL-EHR with built-in Pelstar® wireless technology
1100KLHR	1100KL with mechanical height rod
1100KGHR	1100KG with optional height rod
245EHR-1100	Optional digital height rod for 1100 series*
PROPLUSROD	Optional mechanical height rod for 1100 series
SS-1100KL	ScaleSurance® 2 year extended warranty for 1100 series**
ADPT30	Adapter included, item # for replacement or spare
PLUGSET30	Plug adapters for use outside North America
BTPRINTERKIT-01	Wireless thermal printer kit

Scale also includes connectivity capabilities. To learn about connectivity options, visit [homscales.com/ConnectivitySolutions](https://homscales.com/ConnectivitySolutions)

To learn more, call **1-800-815-6615**  
or visit **[homscales.com/1100KL](https://homscales.com/1100KL)**

Optional mechanical and digital height rods available



**Health o meter®**  
Professional



\* Compatible with 1100 series scales including a date code with the last two digits of 21 or higher. An upgrade option is available for models including a date code with the last two digits of 18 or higher. For assistance, contact Customer Service at 800-815-6615.  
\*\* Not available in all countries.

1100KL\_F&B\_20250825