

THE 1100KL

DIGITAL PLATFORM SCALE WITH EXTRA WIDE HANDRAILS

The Certified ADA Compliant 1100KL is the industry leading bariatric scale that supports up to 1000 pounds and provides the optimum in quality, reliability, and ease of use. The extra wide "live" handrails on both sides of the platform allow the patient to hold on for stability while the scale accurately calculates weight, ideal for use with ambulatory patients. Designed to meet the needs of doctor's offices, hospitals, clinics, and acute and long-term care facilities, the scale's keypad features antimicrobial protection, a display head that swivels 180°, a large slip-resistant low-profile platform, and wheels for easy repositioning.

The scale's high-contrast color screen provides a user friendly interface that includes BMI, hold/release, reweigh, recall, zero, KG/LB, tare, pre-tare and more. The exclusive Everlock® function allows the measuring unit (LB or KG) to be permanently locked on the scale. The scale provides easy integration and seamless communication with today's EMR applications, the 1100KL includes a USB interface for connectivity with a computer, EMR software, or other electronic devices via various Health o meter® Professional Connectivity Solutions.

The 1100KL digital platform scale with extra wide handrails was designed to meet the US Access Boards Standards: Medical facilities operated by state and local governments, also known as public entities, must comply with the MDE Standards for Accessible Medical Diagnostic Equipment, which is enforced by the Department of Justice.

- Slip-resistant platform
- Handrail supports on both sides of platform
- 34" – 38" max. handrail height from platform

Reviewed for ADA compliance and
Platinum certified by Corada

Health o meter®
Professional

Health o meter® Professional 1100KL

Industry leading bariatric scale as reported by Global Industry Analysts, Inc. Patient Scales market survey and Global Healthcare Exchange LLC, (GHX) sales data.



THERE'S STRENGTH IN OUR NUMBERS.

When you add up the numbers on our scale, it's not surprising that the 1100KL is the industry leading bariatric scale.

**100% CERTIFIED
FOR ADA COMPLIANCE**

**120+ LOCAL
FIELD SPECIALISTS**

**100%
CONNECTIVITY**

**#1 TRUSTED
BRAND**

**100 YEARS
OF SERVICE**



weigh**easier**®

TECHNICAL SPECIFICATIONS

Capacity

1000 lb / 454 kg

Resolution

0.2 lb / 0.1 kg

EMR connectivity

Via USB, optional Pelstar® wireless technology

Display

Digit Size: 3/8" / 10 mm
Screen Size: 2 3/4" x 1 1/2" High-Contrast
Color TFT-LCD Screen

Power Source

100-240V adapter (ADPT30)
included, 6 D batteries not included

Functions

LB/KG conversion
LB/KG lock out
BMI
Zero
Tare
Pre-tare
Hold/release
Reweigh
Recall
Auto zero
Auto off
Time/date
Variable auto off time
Audible/mute sound option
Parent/child weighing function
Everlock® - permanent unit lock

Platform Size (w x d x h)

22" x 15 3/4" x 2 1/2"
(559 mm x 400 mm x 64 mm)

Product Footprint (w x d x h)

31 3/4" x 33 1/2" x 49 3/4"
(806 mm x 851 mm x 1264 mm)

Product Weight

52 lb / 23 kg

Portability

2 wheels for ease of movement

Warranty

2-year limited manufacturer warranty

ORDERING INFORMATION

MODEL	DESCRIPTION
1100KL	Digital platform scale with extra wide handrails, LB & KG
1100KG	Digital platform scale with extra wide handrails, KG only
1100KL-BT	1100KL with built-in Pelstar® wireless technology
1100KG-BT	1100KG with built-in Pelstar® wireless technology
1100KL-EHR	1100KL with digital height rod
1100KL-BT-EHR	1100KL-EHR with built-in Pelstar® wireless technology
1100KLHR	1100KL with mechanical height rod
1100KGHR	1100KG with optional height rod
245EHR-1100	Optional digital height rod for 1100 series*
PROPLUSROD	Optional mechanical height rod for 1100 series
SS-1100KL	ScaleSurance® 2 year extended warranty for 1100 series**
ADPT30	Adapter included, item # for replacement or spare
PLUGSET30	Plug adapters for use outside North America
BTPRINTERKIT-01	Wireless thermal printer kit

Scale also includes connectivity capabilities. To learn about connectivity options, visit homscscales.com/ConnectivitySolutions

To learn more, call **1-800-815-6615**
or visit **homscscales.com/1100KL**

Optional mechanical and digital height rods available



Health o meter®
Professional



* Compatible with 1100 series scales including a date code with the last two digits of 21 or higher. An upgrade option is available for models including a date code with the last two digits of 18 or higher. For assistance, contact Customer Service at 800-815-6615.
**Not available in all countries.

1100KL_F&B_20251117