

This classic "Bigfoot" scale has a deep platform with non-slip surface and an 8" easy-to-read raised dial-good for patients with vision challenges or those who prefer more foot room.

## **Advantages**

400 lb / 180 kg Capacity Floor scale The high capacity of this floor scale provides more versatile use

## Weighs in LB/KG

Dial display shows weight in both pounds and kilograms providing versatility

## **8**" Raised Dial Display

Large, easy-to-read dial display is raised above the scale base to allow for more footroom on the platform



**160KL** 

2  $\sigma$ Ο S S Φ Ω

Health o meter® Professional Scales 9500 West 55th Street USA McCook, IL 60525-7110 Phone: 1-800-815-6615 Fax: 1-800-654-7330



## **Specifications**

- Capacity: 400 lb /180 kg
- Graduation: 1 lb / 1 kg
- 8" dial (raised 4")
- Platform Size (w x d x h): 11" x 12 ½" x 3" / 279 mm x 318 mm x 76 mm
- Product Footprint (w x d x h):
  - 11 ¾" x 18 ¼" x 7 ¼" / 298 mm x 464 mm x 184 mm
- Product Weight: 9 lb / 4 kg
- 1 year limited warranty

Model	Description
160KL	Mechanical floor scale - kilograms & pounds (2 pack)
160LB	Mechanical floor scale - pounds only (2 pack)
160KG	Mechanical floor scale - kilograms only (2 pack)
160KLS	Mechanical floor scale - kilograms & pounds (single)
160LBS	Mechanical floor scale - pounds only (single)
160KGS	Mechanical floor scale - kilograms only (single)

