Digital In-Bed / Stretcher Scale

 ϖ O S S Φ Ω



This scale can weigh immobile patients up to 400 lb / 181 kg. It has a durable steel structure, fabric sling, and wheels for portability.

Advantages

EMR Connectivity

Scale has the ability to interface with a computer, EMR software or other electronic device via various Health o meter[®] Professional Connectivity Solutions

Value-Added Features

Operating functions (LB/KG Conversion, LB/KG Lock Out, BMI, Zero, Tare, Pre-tare Hold/Release, Reweigh, Recall, Auto Zero, Auto Off) included at no additional cost

Antimicrobial Keypad The keypad features antimicrobial protection

Calculates Body Mass Index (BMI) Scale functions include the ability to calculate the patient's BMI

Hospital Grade Construction Durable steel structure with large sturdy stretcher enables patients to be weighed accurately and comfortably

Everlock[®] (patent pending) Option that allows users to permanently lock the scale in either LB or KG to improve safety



Digital In-Bed / Stretcher Scale



Manufactured by: Pelstar, LLC. 9500 West 55th St. McCook, IL 60525-7110 USA Phone: 1-800-815-6615 Fax: 1-800-654-7330







Specifications

- Capacity: 400 lb / 181 kg
- Resolution: 0.2 lb / 0.1 kg
- Connectivity via USB, optional Pelstar® wireless technology
- Display: 1 1/2" / 38 mm High-Contrast Color TFT-LCD Screen
- Canvas Size (w x d):
 - 70 ½" x 30"

1791 mm X 762 mm

 Product Footprint - at highest extension of lift (w x d x h): 70" x 44 ½" x 77"

1778 mm x 1130 mm x 1956 mm

- Product Weight: 104 lb / 47 kg
- Functions: LB / KG Conversion, LB / KG lock Out, BMI, Zero, Tare, Pre-tare, Hold / Release, Reweigh, Auto Zero, Auto Off, Everlock®
- 4 Wheels for ease of movement, 1 locking wheel for safety
- Power Source: 100-240V adapter (ADPT30) included,
 - 6 D batteries (not included)
- 2 year limited warranty

Model	Description
2000KL	Digital in-bed / stretcher scale, pounds and kilograms
2000KL-BT	2000KL with Pelstar wireless technology
2000KG	Digital in-bed / stretcher scale, kilograms only
2000KG-BT	2000KG with Pelstar wireless technology
SS-2000KL	ScaleSurance 2 year extended warranty for 2000KL, 2000KG, 2000KL-BT & 2000KG-BT*
ADPT30	Adapter included, item # for replacement or spare
PLUGSET30	Plug adapters for use outside North America
PELSTARONE-P	Wireless thermal printer



* Not available in all countries

2000KL_F&B_20210218