# S S



# **Easy to Use**

Smooth sliding headpiece and telescoping rod allows the user to quickly and easily read the height measurement.

### **Durable Construction**

Aluminum height rod withstands frequent use. Rigorous laboratory tests showed that the Health o meter® Professional 201HR-400 lasts 2 1/2 times longer than similar models.

### **Footpiece Positioner**

Attached footpiece serves as a guide for heel placement and ensures accurate measurements.

### ✓ Large Measurement Range

The 2" / 6cm - 90" / 230cm measuring range provides height readings for a wide range of patients.

Model	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
201HR-400	57" x 3.5" x 3.5" / .4	3.6 lb	854896003870

## **Specifications**

- Measuring Range: 2" / 6cm 90" / 230cm
- Graduation: 1/16" / 1mm
- Physical Dimensions:
  - Width: 2 ½" / 57mm (3" / 76mm including footpiece)
  - Depth: <sup>3</sup>/<sub>4</sub>" / 19mm (1 <sup>3</sup>/<sub>4</sub>" / 44mm including spacers)
  - o Height: 54" / 1370mm H (92" / 2300mm with height rod fully extended)
  - o Headpiece Depth When Unfolded: 11" / 280mm
- 2 Year Warranty

Can be attached to 400KL, 402KL, 402LB, 402KG, and 450KL



Health o meter® Professional Scales 9500 West 55th Street USA McCook, IL 60525-7110 Phone: 1-800-815-6615 Fax: 1-800-654-7330