Wall Mounted Telescopic Metal Height Rod





Easy to Use

Smooth sliding headpiece and telescoping rod allows the user to quickly and easily read the height measurement.

Durable Construction

Aluminum height rod withstands frequent use. Rigorous laboratory tests showed that the Health o meter® Professional 201HR lasts 2 ½ times longer than similar models.

Footpiece Positioner

Attached footpiece serves as a guide for heel placement and ensures accurate measurements.

Large Measurement Range

The 2° / 6cm – 90° / 230cm measuring range provides height readings for a wide range of patients.

Wall Mount Design

Installs easily and quickly with simple tools and provided hardware.



Model	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
201HR	57" x 3.5" x 3.5" / .4	3.6 lb	854896003689

Specifications

- Measuring Range: 2" / 6cm 90" / 230cm
- Graduation: 1/16" / 1mm
- Physical Dimensions:
 - Width: 2 1/4" / 57mm (13" / 330mm including footpiece)
 - Depth: ³/₄" / 19mm (1 ³/₄" / 44mm including spacers)
 - Height: 54" / 1370mm H (92" / 2300mm with height rod fully extended)
 - Headpiece Depth When Unfolded: 11" / 280mm
- 2 Year Warranty

Health o meter[®] Professional www.homscales.com

Health o meter® Professional Scales 9500 West 55th Street USA McCook, IL 60525-7110 Phone: 1-800-815-6615 Fax: 1-800-654-7330

F&B_201HR_20180726