



- ✓ **Easy to Use**
Lockable headpiece and clear indicator allow the user to quickly and easily read the height measurement.
- ✓ **Adjustable Measuring Range**
With a range of over 55 inches, users can select a measuring range that suits their practice, whether they care for children or adults.
- ✓ **Wall Mount Design**
Installs easily and quickly with simple tools and provided hardware.
- ✓ **Durable Construction**
High strength polymer height rod resists bending and ensures that the headpiece slides smoothly.
- ✓ **Choice of Units of Measurement**
Double-sided measuring tape enables users to choose centimeters only or feet/inches/centimeters.



Model	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
205HR	61" x 5" x 4.5" / 0.8	5 lb	854896003696

Specifications

- Height Rod Measurement Range: 55.5" (141 cm),
Measuring Tape Range: 0" - 90" / 0 cm - 230 cm
- Graduation: 1/16" / 1mm
- Physical Dimensions:
 - Width: 2.625" / 6.7 cm
 - Depth: (height rod only): 3/4" / 1.9 cm
 - Height: 60" / 152.4 cm
 - Headpiece Depth: 8" / 20 cm
 - Headpiece Width: 3.875" / 9.8 cm
- 2 Year Warranty

