



- ✓ **Digital Height Measurements**
 Users can quickly and easily read the measurement at any height or angle.
- ✓ **Efficient and Practical Headpiece**
 Well-designed headpiece slides smoothly along height rod channel and remains securely in place when not in use.
- ✓ **Adjustable Measuring Range**
 With a range of over 3 feet (100 cm), users can select a measuring range that suits their practice, whether they care for children or adults.
- ✓ **Wall Mount Design**
 Durable height rod channel installs easily and quickly with simple tools and provided hardware.
- ✓ **Choice of Units of Measurement**
 Digital headpiece provides height measurement in centimeters, feet, or inches.



Specifications

- Height Rod Measurement Range: 39 1/3" / 100 cm (within height rod range of 23 5/8" - 78 3/4" / 60 cm - 200 cm)
- Resolution: 0.1" / 1mm
- Physical Dimensions:
 - Width: 2 1/8" / 54 mm
 - Depth: (rod only): 1/2" / 13 mm
 - Height: 45 5/8" / 1160 mm
 - Headpiece Depth: 7 1/8" / 181 mm
 - Headpiece Width: 1 3/4" / 44 mm
- Power Source: 2 AAA batteries (included)
- Display: 1/2" / 13 mm digit height
- 2 Year Warranty

Model	Description
SS-209HR	ScaleSurance 2 year extended warranty