

## THE 2101KL

### DIGITAL PLATFORM SCALE

The all-digital 2101KL is a comprehensive bariatric scale that supports up to 1000 pounds and provides the optimum in quality, reliability, and ease of use. The extra wide "live" handrails allow the patient to hold on for stability while the scale accurately calculates weight.

Designed to meet the needs of doctor's offices, hospitals, clinics, and acute and long-term care facilities, the scale's keypad features antimicrobial protection, a display head that swivels 180°, a large non-skid low-profile platform, and wheels that enable portability.

The scale's high-contrast color screen provides a user-friendly interface that includes BMI, hold/release, reweigh, recall, zero, KG/LB, tare, pre-tare and more. The exclusive Everlock® function allows the measuring unit (LB or KG) to be permanently locked on the scale.

The scale provides easy integration and seamless communication with today's EMR applications, the 2101KL includes a USB interface for connectivity with a computer, EMR software, or other electronic devices via various Health o meter® Professional Connectivity Solutions.



### THERE'S STRENGTH IN OUR NUMBERS.

When you add up the numbers on our scale, it's not surprising that the 2101KL is your choice for a bariatric platform scale.

**120+** LOCAL  
FIELD SPECIALISTS

**100%**  
CONNECTIVITY

**#1** TRUSTED  
BRAND

**100** YEARS  
OF SERVICE

**Health o meter®**  
Professional

Health o meter® Professional 2101KL

weigh**easier**®

## TECHNICAL SPECIFICATIONS

### Capacity

1000 lb / 454 kg

### Resolution

0.2 lb / 0.1 kg

### EMR Connectivity

Via USB, optional Pelstar® wireless technology

### Display

Digit Size: 3/8" / 10 mm

Screen Size: 2 3/4" x 1 1/2" High-Contrast

Color TFT-LCD Screen

### Power Source

100-240V adapter (ADPT30)

included, 6 D batteries not included

### Functions

LB/KG conversion

LB/KG lock out

BMI

Zero

Tare

Pre-tare

Hold/release

Reweigh

Recall

Auto zero

Auto off

Time/date

Variable auto off time

Audible/mute sound option

Parent/child weighing function

Everlock® - permanent weight lock

### Platform Size (w x d x h)

26" x 22" x 2 3/8"

(660 mm x 559 mm x 60 mm)

### Product Footprint

31" x 27" x 50"

(787 mm x 686 mm x 1270 mm)

### Product Weight

66 lb / 30 kg

### Portability

2 wheels for ease of movement

### Warranty

2-year limited manufacturer warranty

## ORDERING INFORMATION

MODEL	DESCRIPTION
2101KL	Digital platform scale, LB & KG
2101KL-BT	2101KL with Pelstar® wireless technology
2101KG	Digital platform scale, kilograms only
2101KG-BT	2101KG with built-in Pelstar® wireless technology
2101KLHR	2101KL with mechanical height rod included
2101KGHR	2101KG with mechanical height rod included
2101KLHR-BT	2101KL with mechanical height rod included and built-in Pelstar® wireless technology
2101KGHR-BT	2101KG with mechanical height rod included and built-in Pelstar® wireless technology
SS-2101KL	ScaleSurance® 2 year extended warranty*
ADPT30	Adapter included, item # for replacement or spare
PLUGSET30	Plug adapters for use outside North America
245EHR-2101	Optional digital height rod for 2101 models**
STROD	Optional mechanical height rod for 2101 models
BTPRINTERKIT-01	Wireless thermal printer kit

Scale also includes connectivity capabilities. To learn about connectivity options, visit [homescapes.com/ConnectivitySolutions](https://homescapes.com/ConnectivitySolutions)

To learn more, call **1-800-815-6615**  
or visit **[homescapes.com/2101KL](https://homescapes.com/2101KL)**

### Digital and Mechanical Height Rod Options Available



**Health o meter®**  
Professional

\* Not available in all countries

\*\* Compatible with 2101 series scales including a date code with the last two digits of 21 or higher. An upgrade option is available for models including a date code with the last two digits of 17 or higher. For assistance, contact Customer Service at 800-815-6615.

