

THE 2101KL

DIGITAL PLATFORM SCALE

The Certified ADA Compliant 2101 KL is a comprehensive bariatric scale that supports up to 1000 pounds and provides the optimum in quality, reliability, and ease of use. The extra wide “live” handrails on both sides of the platform allow the patient to hold on for stability while the scale accurately calculates weight.

Designed to meet the needs of doctor’s offices, hospitals, clinics, and acute and long-term care facilities, the scale’s keypad features antimicrobial protection, a display head that swivels 180°, a large slip-resistant low-profile platform, and wheels that enable portability.

The scale’s high-contrast color screen provides a user-friendly interface that includes BMI, hold/release, reweigh, recall, zero, KG/LB, tare, pre-tare and more. The exclusive Everlock® function allows the measuring unit (LB or KG) to be permanently locked on the scale.

The scale provides easy integration and seamless communication with today’s EMR applications, the 2101KL includes a USB interface for connectivity with a computer, EMR software, or other electronic devices via various Health o meter® Professional Connectivity Solutions.

The 2101KL digital platform scale was designed to meet the US Access Boards Standards: Medical facilities operated by state and local governments, also known as public entities, must comply with the MDE Standards for Accessible Medical Diagnostic Equipment, which is enforced by the Department of Justice.

- Slip-resistant platform
- Handrail supports on both sides of platform
- 34” – 38” max. handrail height from platform

Reviewed for ADA compliance and
Platinum certified by Corada



REVIEWED FOR ADA COMPLIANCE
CERTIFICATION NO. 5849459

Health o meter®
Professional

Health o meter® Professional 2101 KL

THERE'S STRENGTH IN OUR NUMBERS.

When you add up the numbers on our scale, it's not surprising that the 2101 KL is your choice for a bariatric platform scale.

**100% CERTIFIED
FOR ADA COMPLIANCE**

**120+ LOCAL
FIELD SPECIALISTS**

**100%
CONNECTIVITY**

**#1 TRUSTED
BRAND**

**100 YEARS
OF SERVICE**

weigh**easier**®

TECHNICAL SPECIFICATIONS

Capacity

1000 lb / 454 kg

Resolution

0.2 lb / 0.1 kg

EMR Connectivity

Via USB, optional Pelstar® wireless technology

Display

Digit Size: 3/8" / 10 mm
Screen Size: 2 3/4" x 1 1/2" High-Contrast
Color TFT-LCD Screen

Power Source

100-240V adapter (ADPT30)
included, 6 D batteries not included

Functions

LB/KG conversion
LB/KG lock out
BMI
Zero
Tare
Pre-tare
Hold/release
Reweigh
Recall
Auto zero
Auto off
Time/date
Variable auto off time
Audible/mute sound option
Parent/child weighing function
Everlock® - permanent weight lock

Platform Size (w x d x h)

26" x 22" x 2 3/8"
(660 mm x 559 mm x 60 mm)

Product Footprint

31" x 27" x 50"
(787 mm x 686 mm x 1270 mm)

Product Weight

66 lb / 30 kg

Portability

2 wheels for ease of movement

Warranty

2-year limited manufacturer warranty

ORDERING INFORMATION

MODEL	DESCRIPTION
2101KL	Digital platform scale, LB & KG
2101KL-BT	2101KL with Pelstar® wireless technology
2101KG	Digital platform scale, kilograms only
2101KG-BT	2101KG with built-in Pelstar® wireless technology
2101KLHR	2101KL with mechanical height rod included
2101KGHR	2101KG with mechanical height rod included
2101KLHR-BT	2101KL with mechanical height rod included and built-in Pelstar® wireless technology
2101KGHR-BT	2101KG with mechanical height rod included and built-in Pelstar® wireless technology
SS-2101KL	ScaleSurance® 2 year extended warranty*
ADPT30	Adapter included, item # for replacement or spare
PLUGSET30	Plug adapters for use outside North America
245EHR-2101	Optional digital height rod for 2101 models**
STROD	Optional mechanical height rod for 2101 models
BTPRINTERKIT-01	Wireless thermal printer kit

Scale also includes connectivity capabilities. To learn about connectivity options, visit homescapes.com/ConnectivitySolutions

To learn more, call **1-800-815-6615**
or visit **homescapes.com/2101KL**

Digital and Mechanical Height Rod Options Available



Health o meter®
Professional

* Not available in all countries

** Compatible with 2101 series scales including a date code with the last two digits of 21 or higher. An upgrade option is available for models including a date code with the last two digits of 17 or higher. For assistance, contact Customer Service at 800-815-6615.

