

# THE 2101KL

## DIGITAL PLATFORM SCALE

The Certified ADA Compliant 2101 KL is a comprehensive bariatric scale that supports up to 1000 pounds and provides the optimum in quality, reliability, and ease of use. . The scale's robust 4-point handrail design provides support for patients who need extra stability, allowing them to hold on securely while the scale delivers accurate weight measurements.

Designed to meet the needs of doctor's offices, hospitals, clinics, and acute and long-term care facilities, the scale's keypad features antimicrobial protection, a display head that swivels 180°, a large slip-resistant low-profile platform, and wheels that enable portability.

The scale's high-contrast color screen provides a user-friendly interface that includes BMI, hold/release, reweigh, recall, zero, KG/LB, tare, pre-tare and more. The exclusive Everlock® function allows the measuring unit (LB or KG) to be permanently locked on the scale.

The scale provides easy integration and seamless communication with today's EMR applications, the 2101KL includes a USB interface for connectivity with a computer, EMR software, or other electronic devices via various Health o meter® Professional Connectivity Solutions.

The 2101KL digital platform scale was designed to meet the US Access Boards Standards: Medical facilities operated by state and local governments, also known as public entities, must comply with the MDE Standards for Accessible Medical Diagnostic Equipment, which is enforced by the Department of Justice.

- Slip-resistant platform
- Handrail supports on both sides of platform
- 34" – 38" max. handrail height from platform

Reviewed for ADA compliance and  
Platinum certified by Corada



**Health o meter®**  
Professional

Health o meter® Professional 2101 KL

## THERE'S STRENGTH IN OUR NUMBERS.

When you add up the numbers on our scale, it's not surprising that the 2101 KL is your choice for a bariatric platform scale.

**100% CERTIFIED  
FOR ADA COMPLIANCE**

**120+ LOCAL  
FIELD SPECIALISTS**

**100%  
CONNECTIVITY**

**#1 TRUSTED  
BRAND**

**100 YEARS  
OF SERVICE**

weigh**easier**®



## TECHNICAL SPECIFICATIONS

### Capacity

1000 lb / 454 kg

### Resolution

0.2 lb / 0.1 kg

### EMR Connectivity

Via USB, optional Pelstar® wireless technology

### Display

Digit Size: 3/8" / 10 mm  
Screen Size: 2 3/4" x 1 1/2" High-Contrast  
Color TFT-LCD Screen

### Power Source

100-240V adapter (ADPT30)  
included, 6 D batteries not included

### Functions

LB/KG conversion  
LB/KG lock out  
BMI  
Zero  
Tare  
Pre-tare  
Hold/release  
Reweigh  
Recall  
Auto zero  
Auto off  
Time/date  
Variable auto off time  
Audible/mute sound option  
Parent/child weighing function  
Everlock® - permanent weight lock

### Platform Size (w x d x h)

26" x 22" x 2 3/8"  
(660 mm x 559 mm x 60 mm)

### Product Footprint

31" x 27" x 50"  
(787 mm x 686 mm x 1270 mm)

### Product Weight

66 lb / 30 kg

### Portability

2 wheels for ease of movement

### Warranty

2-year limited manufacturer warranty

## ORDERING INFORMATION

| MODEL           | DESCRIPTION                                                                          |
|-----------------|--------------------------------------------------------------------------------------|
| 2101KL          | Digital platform scale, LB & KG                                                      |
| 2101KL-BT       | 2101KL with Pelstar® wireless technology                                             |
| 2101KG          | Digital platform scale, kilograms only                                               |
| 2101KG-BT       | 2101KG with built-in Pelstar® wireless technology                                    |
| 2101KLHR        | 2101KL with mechanical height rod included                                           |
| 2101KGHR        | 2101KG with mechanical height rod included                                           |
| 2101KLHR-BT     | 2101KL with mechanical height rod included and built-in Pelstar® wireless technology |
| 2101KGHR-BT     | 2101KG with mechanical height rod included and built-in Pelstar® wireless technology |
| SS-2101KL       | ScaleSurance® 2 year extended warranty*                                              |
| ADPT30          | Adapter included, item # for replacement or spare                                    |
| PLUGSET30       | Plug adapters for use outside North America                                          |
| 245EHR-2101     | Optional digital height rod for 2101 models**                                        |
| STROD           | Optional mechanical height rod for 2101 models                                       |
| BTPRINTERKIT-01 | Wireless thermal printer kit                                                         |

Scale also includes connectivity capabilities. To learn about connectivity options, visit [homscalcs.com/ConnectivitySolutions](https://homscalcs.com/ConnectivitySolutions)

To learn more, call **1-800-815-6615**  
or visit **[homscalcs.com/2101KL](https://homscalcs.com/2101KL)**

### Digital and Mechanical Height Rod Options Available



**Health o meter®**  
Professional

\* Not available in all countries

\*\* Compatible with 2101 series scales including a date code with the last two digits of 21 or higher. An upgrade option is available for models including a date code with the last two digits of 17 or higher. For assistance, contact Customer Service at 800-815-6615.

