



The 2400KL is the perfect scale for locations where occasional wheelchair-bound patients need to be weighed. The scale folds easily and requires minimal room for storage. The scale includes a handle and wheels for easy portability.

Advantages

EMR Connectivity

Scale interfaces with a computer, EMR software or other electronic device via various Health o meter[®] Professional Connectivity Solutions

Value-Added Features

Operating functions (LB/KG Conversion, LB/KG Lock Out, BMI, Zero, Tare, Pre-tare, Hold/Release, Reweigh, Recall, Auto Zero, Auto Off, Time/Date, Variable Auto Off Time, Audible/Mute Sound Option) included at no additional cost

Calculates Body Mass Index (BMI) Scale functions include the ability to calculate the patient's BMI

Antimicrobial Keypad The keypad features antimicrobial protection

Remote Display and Adjustable Tracks Display can be mounted on wall or placed on table. Wheelchair rails expand and adjust from 25 ³/₄" to 50 ¹/₈" to accommodate all standard and bariatric wheelchairs

Easy Mobility

Incorporated wheels and casters enable smooth and easy portability, reducing strain on staff

Convenient Storage

Track base folds, can easily be stored in a closet or other small space

Everlock[®] (patent pending)

Option that allows users to permanently lock the scale in either LB or KG to improve safety



Health o meter[®]

www.homscales.com





Specifications

- Capacity: 800 lb / 363 kg
- Resolution: 0.2 lb / 0.1 kg
- Platform Size:
- Expands and adjusts from 25 ¾" (654 mm) to 50 ⅛" (1273 mm)
- Rail Size (w x d):
 - 6" x 40" / 152 mm x 1017 mm
- Product Footprint, Operating Position with Ramps Fully Extended: (w x h x d): 51 ¼" x 40 ½" x 3 ½" / 1302 mm x 1029 mm x 89 mm
- Product Footprint, Folded Position: (w x h x d):
 - 6 ¾" x 40 ½" x 30" / 171 mm x 1029 mm x 762 mm
- Product Weight: 40 lb / 18 kg
- Wheels: (2) Wheels, (2) Casters
- Display: 1 1/2" / 38 mm High-Contrast Color TFT-LCD Screen
- Connectivity: USB, optional Pelstar® wireless technology
- Power Source: 100-240V adapter (ADPT30) included, 6 D batteries (not included)
- Functions: LB/KG Conversion, LB/KG Lock Out, BMI, Zero, Tare, Pre-tare, Hold/Release, Reweigh, Recall, Auto Zero, Auto Off, Time/Date,
- Variable Auto Off Time, Audible/Mute Sound Option, Everlock®
- 2 year limited warranty

Model	Description
2400KL	Digital portable wheelchair scale, pounds & kilograms
2400KL-BT	2400KL with Pelstar [®] wireless technology
2400KG	Digital portable wheelchair scale, kilograms only
2400KG-BT	2400KG with Pelstar [®] wireless technology
SS-2400KL	ScaleSurance 2 year extended warranty (2400 models)*
ADPT30	Adapter included, item # for replacement or spare
PLUGSET30	Plug adapters for use outside North America
PELSTARONE-P	Wireless thermal printer

* Not available in all countries

2400KL_F&B_20210218