

## THE 2400KL

# DIGITAL PORTABLE WHEELCHAIR SCALE

The all-digital 2400KL supports up to 800 pounds and provides the optimum in quality, reliability, and ease of use. The scale folds easily for convenient storage and comes equipped with a handle, wheels, and casters, making it easy to transport and reduces strain on caregivers. The remote display provides versatility with the ability to be mounted on a wall or placed on a table. Designed to meet the needs of hospitals, clinics, and acute and long-term care facilities, the scale's keypad features antimicrobial protection, two ramps to roll the patient on and off the scale, and adjustable rails to accommodate wide wheelchairs.

The scale's high-contrast color screen provides a user friendly interface that includes BMI, hold/release, reweigh, recall, zero, KG/LB, tare, pre-tare and more. The exclusive Everlock® function allows the measuring unit (LB or KG) to be permanently locked on the scale. The scale provides easy integration and seamless communication with today's EMR applications, the 2400KL includes a USB interface for connectivity with a computer, EMR software, or other electronic devices via various Health o meter® Professional Connectivity Solutions.



## THERE'S STRENGTH IN OUR NUMBERS.

When you add up the numbers, it's not surprising that the Health o meter® Professional Scales is the industry-leading name in medical scales.

**120+** LOCAL  
FIELD SPECIALISTS

**2** WHEELS FOR  
EASY PORTABILITY

**100%**  
CONNECTIVITY

**#1** TRUSTED  
BRAND

**100** YEARS  
OF SERVICE

**Health o meter®**  
Professional

Health o meter® Professional 2400KL

weigh**easier**®

# TECHNICAL SPECIFICATIONS

## Capacity

800 lb / 363 kg

## Resolution

0.2 lb / 0.1 kg

## EMR connectivity

Via USB, optional Pelstar® wireless technology

## Display

Digit Size: 3/8" / 10 mm

Screen Size: 2 3/4" x 1 1/2" High-Contrast Color TFT-LCD Screen

## Power Source

100-240V adapter (ADPT30) included, 6 D batteries not included

## Functions

LB/KG conversion  
LB/KG lock out  
BMI  
Zero  
Tare  
Pre-tare  
Hold/release  
Reweigh  
Recall  
Auto zero  
Auto off  
Time/date  
Variable auto off time  
Audible/mute sound option  
Parent/child weighing function  
Everlock®

## Platform Size

Expands and Adjusts from 25 3/4" – 50 1/8"  
(654 mm – 1273 mm)

## Rail Size (w x d)

6" x 40"  
(152 mm x 1017 mm)

## Product Footprint Folded

6 3/4" x 40 1/2" x 30"  
(171 mm x 1029 mm x 762 mm)

## Product Footprint Ramps Fully Extended

51 1/4" x 40 1/2" x 3 1/2"  
(1302 mm x 1029 mm x 89 mm)

## Product Weight

40 lb / 18 kg

## Portability

2 wheels (2 casters)  
for ease of movement

## Warranty

2-year limited manufacturer warranty

## ORDERING INFORMATION

MODEL	DESCRIPTION
2400KL	Digital portable wheelchair scale, LB & KG
2400KG	Digital portable wheelchair scale, KG only
2400KL-BT	2400KL with built-in Pelstar® wireless technology
2400KG-BT	2400KG with built-in Pelstar® wireless technology
SS-2400KL	ScaleSurance® 2 year extended warranty for 2400 series*
ADPT30	Adapter included, item # for replacement or spare
PLUGSET30	Plug adapters for use outside North America
BTPRINTERKIT-01	Wireless thermal printer kit

Scale also includes connectivity capabilities. To learn about connectivity options, visit [hom scales.com/ConnectivitySolutions](http://hom scales.com/ConnectivitySolutions)

To learn more, call 1-800-815-6615  
or visit [hom scales.com/2400KL](http://hom scales.com/2400KL)



**Health o meter®**  
Professional



\* Not available in all countries.

2400KL\_F&B\_20260216