THE 3001 HR

MECHANICAL HEIGHT ROD FOR THE 3001 SERIES OF SCALES

The Health o meter® Professional 3001 HR is an optional mechanical height rod for the Antimicrobial 3001 series of scales. The durable aluminum height rod withstands frequent use, and rigorous laboratory testing shows that the 3001 HR lasts 2 $\frac{1}{2}$ times longer than similar products.

The easy-to-use smooth sliding headpiece and telescoping rod allow the caregiver to read the height measurement quickly and easily. The folding headpiece limits the height rod's profile, providing safety to patients and caregivers. The 3001 HR has an extensive measuring range of 2"-90" in 1/16" increments (6 cm-230cm in 1 mm increments). The optional mechanical height rod has upgraded brackets making assembly quick and easy.



Measuring Range

2" - 90" / 6cm - 230cm Graduation: 1/16" / 1mm

Warranty

2-year limited manufacturer warranty

3001 HR Compatible Scales 3001 KL-AM 3001 KG-AM 3001 KL-AM-BT 3001 KG-AM-BT

ORDERING INFORMATION

MODEL	DESCRIPTION
3001 HR	Optional mechanical height rod for 3001 Series

Health o meter®

To learn more, call 1-800-815-6615 or visit homscales.com/3001HR

weigheasier®