

THE 500HB

HANDLEBAR ACCESSORY FOR THE 500 SERIES OF SCALES

The Health o meter® Professional 500HB is a unique handlebar accessory for the 500KL and is easy to attach to any new scale using the pre-drilled mounting holes and can also be attached to any existing scale using the easy-to-use drill guides. It provides the stability that many patients need without affecting the accuracy of the weight reading.



TECHNICAL SPECIFICATIONS

Product Footprint with 500KL

14 3/8" x 21 1/8" x 53 5/8"
(359 mm x 537 mm x 1362 mm)

Product Weight with 500KL

31 lb / 15 kg

Warranty

2-year limited manufacturer warranty

500HB Compatible Scales

500KL
500KG
500KLAD
500KGAD
500KL-BT
500KG-BT
501KL
501KG
501KL-BT
501KG-BT
502KL
502KG
502KL-BT
502KG-BT

ORDERING INFORMATION

MODEL	DESCRIPTION
500HB	Optional handlebar accessory for 500 series of scales

Health o meter®
Professional

Health o meter® Professional 500HB

weigh**easier**®

To learn more, call **1-800-815-6615**
or visit **homscales.com/500HB**