

## THE 2500KL

# DIGITAL WHEELCHAIR RAMP SCALE

## THE 2500CKL

# DIGITAL WHEELCHAIR DUAL RAMP SCALE WITH FOLD AWAY SEAT

The all-digital 2500KL supports up to 1000 pounds and provide the optimum in quality, reliability, and ease of use. The wide and low profile platform accommodates wide wheelchairs and provides a comfortable weighing experience for patients. Designed to meet the needs of doctor's offices, hospitals, clinics, and acute and long-term care facilities, the scale's keypad features antimicrobial protection, a display head that swivels 180°, a large non-skid low-profile platform, and wheels for easy mobility and allow for convenient placement.

The 2500CKL is perfect for use with standard-size wheelchairs or to safely weigh patients who need to be seated.

The scale's high-contrast color screen provides a user friendly interface that includes BMI, hold/release, reweigh, recall, zero, KG/LB, tare, pre-tare and more. The exclusive Everlock® function allows the measuring unit (LB or KG) to be permanently locked on the scale. The scale provides easy integration and seamless communication with today's EMR applications, the 2500KL and 2500CKL include a USB interface for connectivity with a computer, EMR software, or other electronic devices via various Health o meter® Professional Connectivity Solutions.



2500KL

## THERE'S STRENGTH IN OUR NUMBERS.

When you add up the numbers, it's not surprising that the Health o meter® Professional Scales is the industry-leading name in medical scales.

120+ LOCAL  
FIELD SPECIALISTS

100%  
CONNECTIVITY

#1 TRUSTED  
BRAND

100 YEARS  
OF SERVICE

**Health o meter®**  
Professional

Health o meter® Professional 2500KL & 2500CKL

weigh**easier**®

# TECHNICAL SPECIFICATIONS

## Capacity

1000 lb / 454 kg

## 2500CKL Seat Capacity

419 lb / 190 kg

## Resolution

0.2 lb / 0.1 kg

## EMR connectivity

Via USB, optional Pelstar® wireless technology

## Display

Digit Size: 3/8" / 10 mm  
Screen Size: 2 3/4" x 1 1/2" High-Contrast  
Color TFT-LCD Screen

## Power Source

100-240V adapter (ADPT30) included, 6 D batteries not included

## Functions

LB/KG conversion  
LB/KG lock out  
BMI  
Zero  
Tare  
Pre-tare  
Hold/release  
Reweigh  
Recall  
Auto zero  
Auto off  
Time/date  
Variable auto off time  
Audible/mute sound option  
Parent/child weighing function  
Everlock®

## Platform Size (w x d x h)

29 3/8" x 28 1/4" x 2 3/8"  
(746 mm x 718 mm x 60 mm)

## Ramp Size

28 1/4" x 9 7/8"  
(714 mm x 248 mm)

## 2500CKL Seat Size

18" x 15" x 18 1/2"  
(457 mm x 381 mm x 470 mm)

## Product Weight

2500KL: 81 lb / 37 kg  
2500CKL: 121 lb / 55 kg

## Portability

2 wheels for ease of movement

## Warranty

2-year limited manufacturer warranty

# ORDERING INFORMATION

MODEL	DESCRIPTION
2500KL	Digital wheelchair ramp scale, LB & KG
2500KG	Digital wheelchair dual ramp scale, KG only
2500CKL	Digital wheelchair ramp scale with foldaway seat, LB & KG
2500CKG	Digital wheelchair ramp scale with foldaway seat, KG only
2500KL-BT	2500KL with built-in Pelstar® wireless technology
2500KG-BT	2500KG with built-in Pelstar® wireless technology
SS-2500KL	ScaleSurance® 2 year extended warranty for 2500 series*
SS-2500CKL	ScaleSurance® 2 year extended warranty for 2500CKL & 2500CKG*
CHAIRACC	Optional foldaway chair attachment
B2500RAMP	Optional secondary ramp for 2500 series
ADPT30	Adapter included, item # for replacement or spare
PLUGSET30	Plug adapters for use outside North America
BTPRINTERKIT-01	Wireless thermal printer kit

Scale also includes connectivity capabilities. To learn about connectivity options, visit [homscases.com/ConnectivitySolutions](http://homscases.com/ConnectivitySolutions)

To learn more, call **1-800-815-6615**  
or visit [homscases.com/2500KL](http://homscases.com/2500KL)  
[homscases.com/2500CKL](http://homscases.com/2500CKL)



2500CKL

**Health o meter®**  
Professional



\*Not available in all countries.

\*\*2500CKL does not have MET certification.