

Advantages

EMR Connectivity

Scale has the ability to interface with a computer, EMR software or other electronic device via various Health o meter® Professional Connectivity Solutions

Moveable Arm & Foot Rest

The flip-arm and foot rest allows for easier patient access and patient comfort

Motion Detection

Scale has the ability to weigh accurately despite movement of the subject on the scale

✓ Value-Added Features

Operating functions (LB/KG Conversion, LB/KG Lock Out, Zero Out/Tare, Auto Off, Reweigh) included at no additional cost



9 S S 0



Model	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
2595KL	32.625" x 24.875" x 41.875" / 19.67	100.25 lb	892076002626







Specifications

- Capacity: 600 lb / 272 kg
- Resolution: 0.2 lb / 0.1 kg
- Seat Size: 18 1/4"(w) x 15 3/4" (d) x / 464 mm (w) x 400 mm (d)
- Display: 1 1/8" LCD
- Wheels: 2
- Connectivity: USB
- Power Source: 110-240V Adapter (included),
 6 C-cell batteries (not included)
- Functions: LB / KG Conversion, LB / KG Lock Out,
 Zero Out /Tare, Auto Off



2595KL Display



Health o meter® Professional Scales 9500 West 55th Street McCook, IL 60525-7110 USA Phone: 1-800-815-6615 Fax: 1-800-654-7330