THE 2600KL

DIGITAL WHEELCHAIR SINGLE RAMP SCALE

THE 2610KL

DIGITAL WHEELCHAIR DUAL RAMP SCALE

The all-digital 2600KL is our most popular wheelchair scale that supports up to 1000 pounds and provides the optimum in quality, reliability, and ease of use. The platform folds for convenient storage and the casters (the only scale in the market that has them!) provide easy mobility. The wide and low profile platform accommodates wheelchairs and provides

a comfortable weighing experience for patients. For easy on/ off access the 2610KL offers dual ramps. Designed to meet the needs of doctor's offices, hospitals, clinics, and acute and long-term care facilities, the scale's keypad features antimicrobial protection, a display head that swivels 180°, a large non-skid low-profile platform, and wheels that enable portability.

The scale's high-contrast color screen provides a user friendly interface that includes BMI, hold/ release, reweigh, recall, zero, KG/LB, tare, pretare and more. The exclusive Everlock® function allows the measuring unit (LB or KG) to be permanently locked on the scale. The scale provides easy integration and seamless communication with today's EMR applications, the 2600KL and 2610KL include a USB interface for connectivity with a computer, EMR software, or other electronic devices via various Health o meter® Professional Connectivity Solutions.



THERE'S STRENGTH IN OUR NUMBERS.

When you add up the numbers on our scale, it's not surprising that the 2600KL is our most popular wheelchair scale.

120+ LOCAL FIELD SPECIALISTS

100% Connectivity

#1 TRUSTED BRAND

100 YEARS OF SERVICE

weigheasier

2610KL



Health o meter® Professional 2600KL & 2610KL

TECHNICAL SPECIFICATIONS

Capacity 1000 lb / 454 kg

Resolution 0.2 lb / 0.1 kg

EMR connectivity Via USB, optional Pelstar® wireless technology

Display 1 1⁄2" / 38 mm High-Contrast Color TFT-LCD Screen

Power Source 100-240V adapter (ADPT30) included, 6 D batteries not included

Functions

LB/KG conversion LB/KG lock out BMI Zero Tare Pre-tare Hold/release Reweigh Recall Auto zero Auto off Time/date Variable auto off time Audible/mute sound option Parent/child weighing function Everlock[®]

Platform Size (w x d x h) 36" x 32 ¹/₄" x 2 ¹/₂" (914 mm x 819 mm x 64 mm)

Ramp Size

32 ¹/₄" x 9 ⁷/₈" (819 mm x 248 mm)

2600KL Product Footprint

48 %" x 43 %" x 49 %" (1222 mm x 1114 mm x 1254 mm)

2610KL Product Footprint

56 ¾" x 43 ‰" x 49 %" (1441 mm x 1114 mm x 1254 mm)

Product Weight 2600KL: 93 lb / 42 kg, 2610KL: 98 lb / 44kg

Portability

4 wheels (2 casters) for ease of movement

Warranty

2-year limited manufacturer warranty

ORDERING INFORMATION

MODEL	DESCRIPTION
2600KL	Digital wheelchair single ramp scale, LB & KG
2610KL	Digital wheelchair dual ramp scale, LB & KG
2600KG	Digital wheelchair single ramp scale, KG only
2610KG	Digital wheelchair dual ramp scale, KG only
2600KL-BT	2600KL with built-in Pelstar® wireless technology
2610KL-BT	2610KL with built-in Pelstar® wireless technology
2600KG-BT	2600KG with built-in Pelstar® wireless technology
2610KG-BT	2610KG with built-in Pelstar® wireless technology
B2600RAMP	Optional secondary ramp
SS-2600KL	ScaleSurance® 2 year extended warranty for 2600 series*
SS-2610KL	ScaleSurance® 2 year extended warranty for 2610 series*
ADPT30	Adapter included, item # for replacement or spare
PLUGSET30	Plug adapters for use outside North America
PELSTARONE-P	Wireless thermal printer

To learn more, call **1-800-815-6615** or visit **homscales.com/2600KL homscales.com/2610KL**







Health o meter[®] Professional

