





With its large platform and stainless steel base, the 2842KL can serve in multiple capacities. It can be used as a basic stand-on scale with an extra-large platform, and it can also be used in shipping/receiving or in laundry areas.

Advantages

- ✓ **Durable Construction**
The scale has an easy-to-clean stainless steel solid construction and durable, non-skid mat and can withstand frequent
- ✓ **600 lb Capacity Floor Scale**
The high capacity of this floor scale provides more versatile use
- ✓ **AC Adapter Included**
Scale is powered by a 120V AC adapter ADPT31 (included) or 6 C-Cell batteries (not included) which allows for the convenience of dual power sources
- ✓ **Motion Detection**
Scale has the ability to weigh accurately despite movement of the subject on the scale



Specifications

- Capacity: 600 lb / 270 kg
- Resolution: 0.2 lb / 0.1 kg
- 7/8" / 22 mm LCD display
- Functions: LB / KG Conversion, LB / KG Lock Out, Zero Out, Tare, Auto Off, Variable auto-off Time, Enabling unit view with unit of measure locked, Hold/Release. Reweigh, Pre-Tare, Recall, Everlock®
- Platform Size (w x d x h):
22 1/4" x 42" x 2 1/2" /
565 mm x 1067 mm x 64 mm
- Product Footprint (w x d x h):
22 1/4" x 42" x 2 1/2" /
565 mm x 1067 mm x 64 mm
- Product Weight: 53 lb / 24 kg
- Adapter included (ADPT31), 6 C-cell batteries (not included)
- 2 year limited warranty

Model	Description
SS-2842KL	ScaleSurance 2 year extended warranty for 500 models*
ADPT31	Adapter included, item # for replacement or spare