

Antimicrobial Digital Platform Scale with Extended Handrails

Patented* Antimicrobial Powder Coating uses sustainable silver ions to inhibit the growth and transmission of harmful microbes on treated surfaces.

*Patents US6093407(A), US643241 (B1)

FALL PREVENTION

Safe patient handling and reducing costs associated with inpatient falls are primary concerns for healthcare providers

🧹 "Live" Handrails

Advanced technology allows the patient to hold onto handrails for stability while the scale accurately calculates weight

< Extended Handrails

Longer length handrails provide extra support and help patients maintain balance when stepping onto the platform

Low Profile Platform requires only minimal step up, offering the patient a more comfortable and stable weighing experience

WORKFLOW EFFICIENCY

Simple to use and intuitive scales allow caregivers to focus on patients instead of complicated equipment

🥜 On-Screen Help Menu

Guides users through operating functions and provides easy access to multiple settings

< Interactive Display

High contrast color screen provides a user-friendly interface

EMR Capable

Scale has the ability to interface with a PC or other electronic device via USB or optional Pelstar wireless technology

Everlock[®] (patent pending)

Option that allows users to permanently lock the scale in either LB or KG to improve safety



LIVE HANDRAILS FOR PATIENT STABILITY

INFECTION PREVENTION

Reducing patient's risk against Healthcare Associated Infections is a major challenge for medical facilities

- Antimicrobial* Powder Coating The silver ionic technology suppresses the growth of bacteria and helps protect against cross-contamination
- Antimicrobial Keypad The keypad features antimicrobial protection

DURABILITY

Facilities require a product to withstand frequent use, movement, and disinfection

- Convenient Transport Handle allows the user to easily move the scale to where it is needed in the facility
- Rugged Scuff-Proof Wheels enable smooth and easy portability







(EMR)

Health o meter® Professional Scales 9500 West 55th Street USA McCook, IL 60525-7110 Phone: 1-800-815-6615 Fax: 1-800-654-7330

CONVENIENT TRANSPORT HANDLE





SPECIFICATIONS

- Capacity: 1,000 lb / 454 kg
- Resolution: 0.2 lb / 0.1 kg
- Platform Size (w x h x d):
 22" x 15 ³/₄" x 2 ⁵/₈" / 559 mm x 400 mm x 67 mm
- Handrail Width: 20" / 508 mm
- Product Footprint (w x h x d):
 22" x 34 ½" X 52" / 559 mm x 876 mm x 1321 mm
- Product Weight: 43 lb / 19.5 kg
- Display: 1 1/2" / 38 mm High Contrast Color TFT-LCD Screen
- Rugged Scuff-Proof Wheels: 2
- Power Source: 100-240V Adapter ADPT30 (included), 6 D batteries (not included)
- Connectivity: USB, optional Pelstar wireless technology
- Functions: LB/KG Conversion, LB/KG Lock Out, Body Mass Index (BMI), Zero, Tare, Pre-tare, Hold/Release, Recall, Reweigh, Auto Zero, Auto Off, Time/Date, Variable Auto Off Time, Audible/Mute Sound Option, Everlock[®]
- 2 year limited warranty

Model	Description
3001KL-AMX	Digital platform scale with extended handrails, pounds & kilograms
3001KG-AMX	Digital platform scale with extended handrails, kilograms only
ADPT30	Adapter included, item # for replacement or spare
PLUGSET30	Plug adapters for use outside North America
SS-3001KL-AMX	ScaleSurance 2 year extended warranty for 3001KL-AMX & 3001KG-AMX models*
PELSTARONE-P	Wireless thermal printer



* Not available in all countries

3001KLAMX_E_F&B_20210311