

10 Year Limited
Warranty



Advantages

- ✓ **Steel Frame & Powder Coated Fabrication**
The classic physician scale is built with a robust steel frame and powder-coated column to prevent rust
- ✓ **Exclusive 10 Year Limited Warranty**
Health o meter® Professional's industry leading warranty protects you against defects of materials or workmanship
- ✓ **Height Rod Included**
Includes a convenient height rod so clinicians can measure height and weight all in one place. The easy-to-use height rod provides height measurements from 11 7/8" - 84" / 30cm - 213cm



Model	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
402KL	58.875" x 20.25" x 6.375" / 4.40	33.25 lb	892076002015

Optional Wheels Item # 55000



Optional Counterweights Item # 55090



Specifications

- Capacity: 390 lb / 180 kg
- Graduation: 1/4 lb / 100 g
- Platform Size (w x d x h): 10 1/2" x 14" x 3 1/4" / 267 mm x 356 mm x 83 mm, platform mat: non-slip black
- Fixed Beam, pounds and kilograms can be read on both sides
- Professional Dual Reading Height Rod
Measuring Range : 11 7/8" - 84" / 30cm - 213cm
Graduation: 1/8" / 1 mm
- Optional wheels available (item # 55000)
- Optional counter weights available, to increase capacity to 490lb / 210kg (item # 55090)
- 10 Year Limited Warranty