

10 Year Limited
Warranty



Advantages

- ✓ **Adjustable Unit Measurement Bar**
The mechanical rotating poise bars allow for easy reading of weight in lb or kg
- ✓ **Exclusive 10 Year Limited Warranty**
Health o meter® Professional's industry leading warranty protects you against defects of materials or workmanship
- ✓ **Height Rod Included**
The easy-to-use height rod provides height measurements from 11 $\frac{7}{8}$ " - 84" / 30 cm - 213 cm
- ✓ **Powder Coated Fabrication**
Powder coating inside and out protects it from rust, extending the life of the scale
- ✓ **Engraved Numerals**
Engraved numerals on the poise bars can withstand frequent use, providing long-lasting readability



Exclusive Engraved Numerals

Easy to read rotating bars



Powder coated inside and out for extended life

Easy to read rotating bars



Exclusive rotating poise bar

Model	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
450KL	58.875" x 20.25" x 6.375" / 4.40	36.75 lb	892076002169

Specifications

- Capacity: 500 lb / 200 kg – with included counterweights.
- Graduation: 1/4 lb / 100 g
- Platform Size: 10 1/2" (w) x 14" (d) x 3 1/4" (h) / 267 mm (w) x 356 mm (d) x 83 mm (h)
- Professional Dual Reading Height Rod
 - Measuring Range : 11 7/8" - 84" / 30cm - 213cm
 - Graduation: 1/8" / 1 mm
- Wheels: Optional item #55000
- 10 Year Limited Warranty

Optional Wheels Item # 55000

