



The 499KLHR is a compact version of our most popular scale, the 500KL. The waist-high display is easy to read and the scale is an attractive, robust option for use in medical offices. With a height rod included, the 499KLHR is an economical choice for almost any environment.

Advantages

- EMR Connectivity**
Scale has the ability to interface with a computer, EMR software or other electronic device via various Health o meter® Professional Connectivity Solutions
- Calculates Body Mass Index (BMI)**
Scale functions include the ability to calculate the patient's BMI
- Low Profile Platform**
Low platform height requires only minimal step up, offering patient a more comfortable and stable weighing experience
- Value-Added Features**
Operating functions (LB/KG Conversion, LB/KG Lock Out, BMI, Zero, Tare, Hold/Release, Auto Zero, Auto Off) included at no additional cost
- Compact Size**
The slim design of the scale and waist-high display require a limited amount of space



Specifications

- Capacity: 500 lb / 220 kg
- Resolution: 0.2 lb / 0.1 kg
- EMR connectivity via USB, optional Pelstar® wireless technology
- 1" / 25 mm, LCD display
- Functions: LB / KG Conversion, LB / KG Lock Out, BMI, Zero, Tare, Hold / Release, Auto Zero, Auto Off
- Platform Size (w x d x h):
13 3/4" x 16 1/2" x 2 3/8" / 349 mm x 419 mm x 60 mm
- Product Footprint (w x d x h):
14 1/8" x 21 1/8" x 49" / 359 mm x 537 mm x 1245 mm
- Product weight: 22.4 lb / 10 kg
- Height rod: 11 7/8" - 84" / 30 cm - 213 cm
- Height rod graduation: 1/8" / 1 mm
- 2 wheels for ease of movement
- 6 AA batteries included, optional adapter (order # ADPT31)
- 2 year limited warranty

| Model | Description |
|--------------|-------------------------------------------------------------------|
| 499KLHR | Waist high digital scale, pounds & kilograms, height rod included |
| 499KGHR | Waist high digital scale, kilograms only, height rod included |
| 499KLHRAD | 499KL with height rod & adapter ADPT31 included |
| 499KGHRAD | 499KG with height rod & adapter ADPT31 included |
| SS-499KL | ScaleSurance 2 year extended warranty for 499 scale models* |
| ADPT31 | Optional adapter |
| PLUGSET30 | Plug adapters for use outside North America |
| PELSTARONE-P | Wireless thermal printer |