

## THE 499KLHR

# WAIST-HIGH PRIMARY CARE DIGITAL SCALE WITH HEIGHT ROD

The 499KLHR is a compact version of our number-one-selling physician scale in North America the 500KL primary care scale. The scale provides the optimum in quality, reliability, and ease of use. Designed to meet the needs of doctor's offices, hospitals, clinics, acute and long-term care facilities, scale features a slim compact design, a large non-skid platform, and wheels that enable portability. The 499KLHR gives providers and patients a waist-high weight indicator for patients up to 460 pounds (209 kilograms), with a user-friendly digital interface that includes BMI, hold/release, zero, tare, and KG/LB conversion. The exclusive Everlock® function allows the measuring unit (LB or KG) to be permanently locked on the scale

With its low platform height, the 499KLHR is easy to set up and offers patients a more comfortable and stable weighing experience.

### THERE'S STRENGTH IN OUR NUMBERS.

When you add up the numbers on our scale, it's not surprising that the 499KLHR is reliable and easy to use.

**150,000+**  
INSTALLED  
NATIONWIDE

**120+** LOCAL  
FIELD SPECIALISTS

**8" WIDE WHEELBASE**  
BEST IN ITS CLASS

**100 YEARS**  
OF SERVICE

weigh**easier**®

**Health o meter®**

Professional

Health o meter® Professional 499KLHR



## TECHNICAL SPECIFICATIONS

### Capacity

460 lb / 209 kg

### Resolution

0.2 lb / 0.1 kg

### Display

1" / 25 mm, LCD display

### Batteries

6 AA batteries included,  
optional adapter  
(order# ADPT31)

### Functions

LB/KG conversion  
LB/KG lock out  
BMI  
Zero  
Tare  
Hold/Release  
Auto Zero  
Auto Off  
Parent/Child weighing function  
Everlock® - permanent unit lock

### Platform Size (w x d x h)

13  $\frac{3}{4}$ " x 16  $\frac{1}{2}$ " x 2  $\frac{3}{8}$ "  
(349 mm x 419 mm x 60 mm)  
Dent and rust proof platform

### Product Footprint

14  $\frac{1}{8}$ " x 21  $\frac{1}{8}$ " x 43  $\frac{5}{8}$ "  
(359 mm x 537 mm x 1108 mm)

### Height Rod

11" - 84" / 30 cm - 213 cm  
Height rod graduation:  $\frac{1}{8}$ " / 1 mm

### Product Weight

22 lb / 10 kg

### Portability

2 wheels, wide 8" wheelbase for ease  
of movement for ease of movement

### Warranty

2-year limited manufacturer warranty

## ORDERING INFORMATION

MODEL	DESCRIPTION
499KLHR	499KL with height rod included
499KGHR	499KG with height rod included
SS-499KL	ScaleSurance® 2 year extended warranty for 499 series*
ADPT31	Optional adapter
PLUGSET30	Plug adapters for use outside North America

To learn more, call **1-800-815-6615**  
or visit **hom scales.com/499KLHR**



**Health o meter®**  
Professional



\*Not available in all countries.