

THE 501KL

TAA COMPLIANT PRIMARY CARE DIGITAL SCALE

The TAA Compliant number-one-selling physician scale in North America, the all-digital Health o meter® 501KL primary care scale provides the optimum in quality, reliability, and ease of use. Designed to meet the needs of doctor's offices, hospitals, clinics, and acute and long-term care facilities, the scale features a built-in height rod in the scale column, a large non-skid platform, and wheels that enable portability. The scale gives providers and patients an eye-level weight indicator for patients up to 550 pounds, with a user-friendly digital interface that includes BMI, hold/release, zero, tare, reweigh, recall, KG/LB, clear, and enter functions. The exclusive Everlock® function allows the measuring unit (LB or KG) to be permanently locked on the scale.

The first healthcare scale to provide easy integration and seamless communication with today's EMR applications, the 501KL includes a USB interface for connectivity with a computer, EMR software, or other electronic devices via various Health o meter® Professional Connectivity Solutions.

With its low platform height, the 501KL is easy to set up and offers patients a more comfortable and stable weighing experience.

500HB "live" wrap-around handlebar

The 500HB is a unique handlebar accessory for the 500KL and is easy to attach to any new scale using the pre-drilled mounting holes and can also be attached to any existing scale using the easy-to-use drill guides. It provides the stability that many patients need without affecting the accuracy of the weight reading.

Health o meter®
Professional

Health o meter® Professional 501KL

THERE'S STRENGTH IN OUR NUMBERS.

The 501KL is the newest generation of the 500KL, the most widely used scale in physician offices nationwide.

150,000+
INSTALLED
NATIONWIDE

120+ LOCAL
FIELD SPECIALISTS

100%
CONNECTIVITY

#1 TRUSTED
PHYSICIAN SCALE

100 YEARS
OF SERVICE

weigh**easier**®



TECHNICAL SPECIFICATIONS

Capacity

550 lb / 250 kg

Resolution

0.2 lb / 0.1 kg

EMR connectivity

Via USB, optional Pelstar® wireless technology

Display

1" / 25 mm, LCD display

Batteries

6 AA batteries included, optional 100V-240V adapter (order# ADPT30)

Functions

LB/KG conversion
LB/KG lock out
BMI
Zero
Tare, Pre-tare
Hold/Release
Auto Zero
Auto Off
Reweigh
Recall
Parent/Child weighing function
Everlock®

Platform Size (w x d x h)

13 3/4" x 16 1/2" x 2 3/8"
(349 mm x 419 mm x 60 mm)
Dent and rust proof platform

Product Footprint

14 3/8" x 21 1/8" x 53 5/8"
(359 mm x 537 mm x 1362 mm)

Product Weight

21 lb / 9 kg

Height Rod

30"-84 1/4" / 76 cm-214 cm
Graduation: 1/8" / 1 mm

Portability

2 wheels, wide 8" wheelbase for ease of movement

Warranty

2-year limited manufacturer warranty

ORDERING INFORMATION

MODEL	DESCRIPTION
501KL	TAA Compliant eye-level digital scale, LB & KG
501KG	TAA Compliant eye-level digital scale, KG only
501KLAD	501KL with adapter ADPT30
501KGAD	501KG with adapter ADPT30
501KL-BT	501KL with built-in Pelstar® wireless technology
501KG-BT	501KG with built-in Pelstar® wireless technology
500HB	Optional handlebar accessory for 501 series scales
SS-501KL	ScaleSurance® 2 year extended warranty for 501 series*
ADPT30	Optional 100V-240V adapter
PLUGSET30	Plug adapters for use outside North America
PELSTARONE-P	Thermal wireless printer

Model **501KL** with
Optional Handlebar
Accessory **500HB**



To learn more, call **1-800-815-6615**
or visit **homscapes.com/501KL**

Health o meter®

Professional

COUNTRY OF ORIGIN: USA



*Not available in all countries.

501KL_F&B_20201208