

## THE 501KL

# TAA COMPLIANT PRIMARY CARE DIGITAL SCALE

The TAA Compliant number-one-selling physician scale in North America, the all-digital Health o meter® 501KL primary care scale provides the optimum in quality, reliability, and ease of use. Designed to meet the needs of doctor's offices, hospitals, clinics, and acute and long-term care facilities, the scale features a built-in height rod in the scale column, a large non-skid platform, and wheels that enable portability. The scale gives providers and patients an eye-level weight indicator for patients up to 550 pounds, with a user-friendly digital interface that includes BMI, hold/release, zero, tare, reweigh, recall, KG/LB, clear, and enter functions. The exclusive Everlock® function allows the measuring unit (LB or KG) to be permanently locked on the scale.

The first healthcare scale to provide easy integration and seamless communication with today's EMR applications, the 501KL includes the option to add wireless connectivity to interface with a computer, EMR software, or other electronic device.

With its low platform height, the 501KL is easy to set up and offers patients a more comfortable and stable weighing experience.

### 500HB "live" wrap-around handlebar

The 500HB is a unique handlebar accessory for the 500 series and is easy to attach to any new scale using the pre-drilled mounting holes and can also be attached to any existing scale using the easy-to-use drill guides. It provides the stability that many patients need without affecting the accuracy of the weight reading.

## THERE'S STRENGTH IN OUR NUMBERS.

The 501KL is the newest generation of the 500KL, the most widely used scale in physician offices nationwide.

**150,000+**  
INSTALLED  
NATIONWIDE

**120+** LOCAL  
FIELD SPECIALISTS

**#1** TRUSTED  
PHYSICIAN SCALE


**100 YEARS**  
OF SERVICE

**Health o meter®**

**P r o f e s s i o n a l**

Health o meter® Professional 501KL

weigh**easier**®



## TECHNICAL SPECIFICATIONS

### Capacity

550 lb / 250 kg

### Resolution

0.2 lb / 0.1 kg

### EMR connectivity

Via optional Pelstar®  
wireless technology

### Display

1" / 25 mm, LCD display

### Batteries

6 AA batteries included,  
optional 100V-240V adapter  
(order# ADPT30)

### Functions

LB/KG conversion  
LB/KG lock out  
BMI  
Zero  
Tare, Pre-tare  
Hold/Release  
Auto Zero  
Auto Off  
Reweigh  
Recall  
Parent/Child weighing function  
Everlock®

### Platform Size (w x d x h)

13 3/4" x 16 1/2" x 2 3/8"  
(349 mm x 419 mm x 60 mm)  
Dent and rust proof platform

### Product Footprint

14 3/8" x 21 1/8" x 53 5/8"  
(359 mm x 537 mm x 1362 mm)

### Product Weight

21 lb / 9 kg

### Height Rod

30"-84 1/4" / 76 cm-214 cm  
Graduation: 1/8" / 1 mm

### Portability

2 wheels, wide 8" wheelbase for ease  
of movement

### Warranty

2-year limited manufacturer warranty

## ORDERING INFORMATION

MODEL	DESCRIPTION
501KL	TAA Compliant eye-level digital scale, LB & KG
501KG	TAA Compliant eye-level digital scale, KG only
501KLAD	501KL with adapter ADPT30
501KGAD	501KG with adapter ADPT30
501KL-BT	501KL with built-in Pelstar® wireless technology
501KG-BT	501KG with built-in Pelstar® wireless technology
500HB	Optional handlebar accessory for 501 series scales
SS-501KL	ScaleSurance® 2 year extended warranty for 501 series*
ADPT30	Optional 100V-240V adapter
PLUGSET30	Plug adapters for use outside North America

Model **501KL** with  
Optional Handlebar  
Accessory **500HB**



To learn more, call **1-800-815-6615**  
or visit **homscaler.com/501KL**

**Health o meter®**  
Professional  
**COUNTRY OF ORIGIN: USA**



\*Not available in all countries.