Digital Pediatric Seat Scale



The deep infant seat on the 524KL reduces the movement of young children to speed up the weighing process. Weight is captured easily while the seat supports the child's head and body.

Advantages

EMR Connectivity

 ϖ

S

S

Φ

Scale has the ability to interface with a computer, EMR software or other electronic device via various Health o meter[®] Professional Connectivity Solutions

1" / 25 mm LCD Display

The large display shows the user an easy-to-read weight measurement

Motion Detection

Scale has the ability to weigh accurately despite movement of the subject on the scale

Value-Added Features

Operating functions (LB/KG Conversion, LB/KG Lock Out, Zero Out/Tare, Hold/Release, Auto Zero, Auto Off) included at no additional cost

Convenient Weighing Seat

Plastic, easy-to-clean seat limits infant movement while accurately calculating weight measurement



See reverse side for additional specifications

Digital Pediatric Seat Scale



Health o meter® Professional Scales 9500 West 55th Street USA McCook, IL 60525-7110 Phone: 1-800-815-6615 Fax: 1-800-654-7330



Specifications

EMR

- Capacity: 50 lb / 23 kg
- Resolution: 0 20 lb 0.2 oz / 20 50 lb 0.5 oz 0 9 kg / 5g; 9 - 23 kg / 10 g
- Connectivity via USB, optional Pelstar® wireless technology
- 1" / 25 mm LCD display
- Functions: LB/KG Conversion, LB/KG Lock Out, Zero, Tare, Hold / Release, Auto Zero, Auto Off
- Seat Size (w x d x h): 24 ¼" x 14 ¼" x 18 ½" 616 mm x 362 mm x 470 mm
- Product Footprint (w x d x h): 24 ½" x 21" x 23 ½" 622 mm x 533 mm x 587 mm
- Product Weight: 21 lb / 10 kg
- 6 AA batteries included, Optional adapter (order # ADPT31)
- 2 year limited warranty

Model	Description
524KL	Digital pediatric seat scale, pounds & kilograms
524KG	Digital pediatric seat scale, kilograms only
SS-524KL	ScaleSurance 2 year extended warranty
ADPT31	Optional adapter
PELSTARONE	External wireless kit
PELSTARONE-P	Wireless thermal printer

