



This scale is ideal for infants and with its digital height rod practitioners can measure both weight and height simultaneously.

### Advantages

- ✓ **EMR Connectivity**  
Scale has the ability to interface with a computer, EMR software or other electronic device via various Health o meter® Professional Connectivity Solutions
- ✓ **Scale Value-Added Features**  
Operating functions (LB/KG Conversion, KG Lock Out, Zero Out/Tare, Hold/Release, Auto Off) included at no additional cost
- ✓ **Motion Detection**  
Scale has the ability to weigh accurately despite movement of the subject on the scale
- ✓ **Digital Height Rod Included**  
Easy to read height rod allows user to measure both weight and height simultaneously



### Scale Specifications

- Capacity: 44 lb / 20 kg
- Resolution: 0.5 oz / 10 g
- Tray Size (w x d x h):  
24 7/8" x 13" x 1 3/8"  
632 mm x 330 mm x 35 mm
- Display: 1" / 25 mm LCD
- Power Source: 6 AA batteries included,  
optional power adapter  
(order # ADPT30)
- Connectivity: USB, optional Pelstar®  
wireless technology
- Functions: Zero, Tare, Hold / Release,  
Auto Off
- 2 year limited warranty

### Height Rod Specifications

- Measurement Range:  
13 3/4" – 31 1/2" / 35 -80 cm
- Resolution: 1/16" / 0.1 cm
- Power Source: 2 AAA batteries (included)
- Display: LCD, 5/8" / 16 mm digit height
- Digital display provides height in  
centimeters, inches, & feet
- Functions: Auto Off, Hold
- 2 year limited warranty

Model	Description
553KL-EHR	Digital pediatric scale (pounds & kilograms) with digital height rod
553KG-EHR	Digital pediatric scale (kilograms only) with digital height rod
2210CART	Rolling cart for 553KL-EHR & 553KG-EHR
ADPT30	Optional 100-240V adapter for scale
PLUGSET30	Plug adapters for use outside North America for scale
PELSTARONE-P	Wireless thermal printer for scale