

THE 553KL-EHR

DIGITAL PEDIATRIC SCALE WITH EXTRA-WIDE TRAY AND DIGITAL HEIGHT ROD

The 553KL-EHR pediatric scale was specifically created for nurseries, delivery rooms, and pediatric practices, making it an ideal choice for newborns, infants, and toddlers. With a maximum capacity of 44 pounds and 20 kilograms, this scale offers providers and patients a range of benefits, including safety, comfort, and flexibility. The easy to read digital height rod makes measuring infant height fast and easy with clear digital measurements, large head and foot paddles, and convenient Hold function. The scale ensures a secure platform, keeps patients calm for accurate readings, and accommodates infants lying down to toddlers sitting up.

With its user-friendly digital interface, the 553KL scale includes features such as Zero, Tare, Pre-Tare, Hold/Release, Reweigh, and Auto-off, making it convenient and easy to use. Its motion detection technology ensures accurate weight measurements even when the patient moves.

The 553KL-EHR provides easy integration and seamless communication with today's EMR applications with the option to incorporate wireless connectivity for interfacing with computers, EMR software, or other electronic devices.



THERE'S STRENGTH IN OUR NUMBERS.

When you add up the numbers on our scale, it's not surprising that the 553KL-EHR is your choice for safety, comfort, and reliability.

120+ LOCAL
FIELD SPECIALISTS

30-60 SEC AUTO
WEIGHT HOLD

#1 TRUSTED
BRAND

100 YEARS
OF SERVICE

Health o meter[®]
P r o f e s s i o n a l

Health o meter[®] Professional 553KL-EHR

weigh**easier**[®]

TECHNICAL SPECIFICATIONS

Capacity

44 lb / 20 kg

Resolution

0.5 oz / 10g

EMR Connectivity

Via optional Pelstar® wireless technology

Display

1" / 25 mm, LCD display

Power Source

6 AA batteries included, 100-240V adapter (ADPT31) not included

Functions

LB / KG conversion
LB / KG lock out
Zero
Tare
Pre-Tare
Hold / release
Reweigh
Recall
Auto zero
Auto off
Motion detection
Variable auto-off time
Enabling unit view with unit of measure locked
Auto-weight hold for 30 to 60 seconds
Everlock® - permanent weight lock

Tray Size (w x d x h)

24 7/8" x 13" x 1 3/8"
(632 mm x 330 mm x 35 mm)

Product Footprint

24 7/8" x 16 1/2" x 3 7/8"
(632 mm x 419 mm x 98 mm)

Product Weight

10 lb / 5 kg

Digital Height Range

13 3/4" - 31 1/2" / 35 cm - 80 cm
Resolution: 1/16" / 0.1 cm

Height Rod Power Source

2 AAA batteries (included)

Height Rod Display

LCD, 5/8" / 16 mm digit height

Portability

2210CART- optional rolling cart

Warranty

2-year limited manufacturer warranty

ORDERING INFORMATION

MODEL	DESCRIPTION
553KL-EHR	Digital Pediatric Scale with Extra-Wide Tray and Digital Height Rod, LB & KG
553KG-EHR	Digital Pediatric Scale with Extra-Wide Tray and Digital Height Rod, KG only
SS-553KL	ScaleSurance® 2 year extended warranty for 553 series*
2210CART	Rolling scale/ medical device cart for 553 series of scales
ADPT31	Optional adapter for 553 series of scales
ADPT30WA	Optional adapter and plugset for use outside North America
BTPRINTERKIT-01	Wireless thermal printer kit

Scale also includes connectivity capabilities. To learn about connectivity options, visit homscales.com/ConnectivitySolutions



2210CART
Medical Device Cart for 553
Series of Scales

To learn more, call **1-800-815-6615**
or visit homscales.com/553KL-EHR

Health o meter®
Professional



*Not available in all countries.