# Digital Pediatric Tray Scale with Mechanical Height Rod

Order 553KG-HR for kilograms only scale

553KL-HR



This scale is ideal for infants and with its mechanical height rod practitioners can measure both weight and height simultaneously.

# **Advantages**

#### **EMR** Connectivity

Scale has the ability to interface with a computer, EMR software or other electronic device via various Health o meter<sup>®</sup> Professional Connectivity Solutions

#### Scale Value-Added Features

Operating functions (LB/KG Conversion, KG Lock Out, Zero Out/Tare, Hold/Release, Auto Off) included at no additional cost

#### Motion Detection

Scale has the ability to weigh accurately despite movement of the subject on the scale

#### Mechanical Height Rod Included

Easy to read height rod allows user to measure both weight and height simultaneously



See reverse side for additional specifications

Order 553KG-HR for kilograms only scale

Health o meter® Professional Scales

9500 West 55th Street USA

McCook, IL 60525-7110

Phone: 1-800-815-6615

Fax: 1-800-654-7330



### **Scale Specifications**

(EMR)

- Capacity: 44 lb / 20 kg
- Resolution: 0.5 oz / 10 g
- Tray Size (w x d x h): 24 <sup>7</sup>/<sub>6</sub>" x 13" x 1 <sup>3</sup>/<sub>6</sub>" 632 mm x 330 mm x 35 mm
- Display: 1" / 25 mm LCD
- Power Source: 6 AA batteries included, optional power adapter (order # ADPT30)
- Connectivity: USB, optional Pelstar<sup>®</sup> wireless technology
- Functions: Zero, Tare, Hold / Release, Auto Off
- 2 year limited warranty

## **Height Rod Specifications**

- Range of measurement: 7 cm – 100 cm / 3" – 39"
- Division of measurement: 1 mm / 1/16"

Model	Description
553KL-HR	Digital pediatric scale (pounds & kilograms) with mechanical height rod
553KG-HR	Digital pediatric scale (kilograms only) with mechanical height rod
2210CART	Rolling cart for 553KL-HR & 553KG-HR
ADPT30	Optional 100-240V adapter for scale
PLUGSET30	Plug adapters for use outside North America for scale
PELSTARONE-P	Wireless thermal printer for scale



\* Not available in all countries

553KL-HR\_F&B\_20190828