



Advantages

- ✓ **Easy to Move**
Lightweight frame with 4 casters make portability easy for bringing the scale to the patient's location.
- ✓ **High Capacity**
High 600 lb / 270 kg capacity provides flexibility to weigh almost any patient.
- ✓ **Designed for Patient Comfort**
Large seat with moving armrests allow for easy access on and off the scale.
- ✓ **Convenient Storage Features**
Folding seat back and armrests make the scale simple to store.
- ✓ **Value-Added Features**
Operating functions (LB/KG Conversion, LB/KG Lock Out, Body Surface Area (BSA), BMI, Zero, Tare, Hold/Release, Auto Zero, Auto Off) included at no additional cost



Specifications

- Capacity: 600 lb / 270 kg
- Resolution: 0.2 lb / 0.1 kg
- Seat Size (w x d):
18" x 18 ½"
457 mm x 470 mm
- Product Footprint (w x d x h):
23 ¼" x 40 ⅝" x 37 ⅞"
590 mm x 1020 mm x 960 mm
- Product Weight: 39 lb / 18 kg
- Display: 1 ⅜" / 35 mm LCD
- Wheels: 4 Casters for ease of movement
- Power Source: 6 AA batteries (included)
- Functions: LB/KG Conversion, LB/KG Lock Out, BSA, BMI, Zero, Tare, Hold/Release, Auto Zero, Auto Off,
- 2 year limited warranty

Model	Description
SS-594KL	ScaleSurance 2 year extended warranty