



With a high capacity and mechanical height rod, this scale properly weighs and measures the heaviest ambulatory patients.

The 350° swivel head makes it easy for clinicians to view measurements and interact with the scale.

Advantages

EMR Capable

Scale interfaces with a computer, EMR software or other electronic device via Health o meter® Professional Connectivity Solutions

350° Swivel Display 350° Swivel Display

Eye-level rotating display can be read from almost any angle and allows for ease of use and patient privacy

▼ Built-in Height Rod

The easy-to-use height rod provides height measurements from 11 \(\frac{1}{8} \) - 84" / 30 cm - 213 cm

▼ Value-Added Features

Operating functions (LB/KG Conversion, LB/KG Lock Out, BMI, Zero, Tare, Hold/Release, Auto Zero, Auto Off) included at no additional cost

Calculates Body Mass Index (BMI)
Scale functions include the ability to

Scale functions include the ability to calculate the patient's BMI



$\boldsymbol{\omega}$ S S (1)







Specifications



Capacity: 660 lb / 300 kgResolution: 0.2 lb / 0.1 kg

Platform Size (w x h x d):
 14 ¼" x 14 ¼" x 2 5%"
 362 mm x 362 mm x 67 mm

Product Footprint (w x h x d):
 14 ¾" x 23 ¼" x 54 ½"
 375 mm x 591 mm x 1384 mm

• Product Weight: 24 lb / 11 kg

Height Rod: 11 ½" - 84" / 30 cm - 213 cm
Display: 1" / 25 mm LCD, 350° Swivel

• Wheels: 2

Connectivity: USB, optional Pelstar wireless technology

Power Source: 100-240V adapter (ADPT50) included,
 6 C-cell batteries (not included)

Functions: LB / KG Conversion, LB / KG Lock Out , BMI, Zero,
 Tare, Hold / Release, Auto Zero, Auto Off

• 2 year limited warranty

Model	Description
597KG	Heavy duty eye level digital scale, kilograms only
SS-597KL	ScaleSurance 2 year extended warranty for 597 models*
ADPT50	Adapter included, item # for replacement or spare
PLUGSET50	Plug adapters for use outside North America
PELSTARONE-P	Wireless thermal printer

Health o meter® Professional Scales 9500 West 55th Street USA McCook, IL 60525-7110 Phone: 1-800-815-6615

Fax: 1-800-654-7330

* Not available in all countries

