Heavy Duty Eye Level Digital Scale with Integral Digital Height Rod



Improve workflow with this top-of-the-line physician office scale. The built-in digital height rod eliminates the need to manually enter height for computing BMI. The 350° swivel head allows for line-of-sight viewing and patient privacy.

Advantages

EMR Capable

Scale interfaces with a computer, EMR software or other electronic device via Health o meter[®] Professional Connectivity Solutions

350° Swivel Display

Eye-level rotating display can be read from almost any angle and allows for ease of use and patient privacy

Built-in Digital Height Rod

The easy-to-use, state of the art digital height rod provides measurements from $23 \frac{5}{8}$ " - 88" / 60 cm - 222 cm

Value-Added Features

Operating functions (LB/KG Conversion, LB/KG Lock Out, BMI, Zero, Tare, Hold/Release, Auto Zero, Auto Off) included at no additional cost

Calculates Body Mass Index (BMI)

Scale functions include the ability to calculate the patient's BMI

Patented Design

The 600KL's patented design provides a user-friendly and intuitive interface for measuring patient weight, height and BMI with a convenient rotating display head

See reverse side for additional specifications



Heavy Duty Eye Level Digital Scale with Integral Digital Height Rod



Specifications

- Capacity: 600 lb / 272 kg
- Resolution: 0.2 lb / 0.1 kg
- Height Rod Resolution: 0.1" / 0.1cm
- Platform Size (w x d x h): 14 ¼" x 14 ¼" x 2 5%" 362 mm x 362 mm x 67 mm

US Patent #D667323

EMF

NE.

- Product Footprint (w x h x d): 14 ³⁄₄" x 24 ¹⁄₆" x 54 ³⁄₈" 375 mm x 613 mm x 1381 mm
- Product Weight: 28 lb / 13 kg
- Digital Height Rod: 23 5/8" 88" / 60 cm 222 cm
- Connectivity: USB, optional Pelstar wireless technology
- 1" / 25 mm LCD display, 35° left / right swivel (350° with height rod fully retracted) for ease of viewing
- Functions: LB / KG Conversion, LB / KG Lock Out, BMI, Zero, Tare, Hold / Release, Auto Zero, Auto Off
- 2 wheels for ease of movement
- 100 240V adapter (ADPT50) included, 6 C-cell batteries (not included)
- 2 year limited warranty

Model	Description
600KG	Heavy duty eye level digital scale with integral digital height rod, KG only
600KL-BT	600KL with built-in Pelstar wireless technology
600KG-BT	600KG with built-in Pelstar wireless technology
SS-600KL	ScaleSurance 2 year extended warranty for 600KL & 600KG*
SS-600KL-BT	ScaleSurance 2 year extended warranty for 600KL-BT & 600KG-BT*
ADPT50	Adapter included, item # for replacement or spare
PLUGSET50	Plug adapters for use outside North America
PELSTARONE-P	Wireless thermal printer

Health o meter® Professional Scales 9500 West 55th Street USA McCook, IL 60525-7110 Phone: 1-800-815-6615 Fax: 1-800-654-7330

(D)


Ο

S

S

Φ

Ω



600KL_F&B_20210412

Health o meter® Professional www.homscales.com