



US Patent #D667323

Improve workflow with this top-of-the-line physician office scale. The built-in digital height rod eliminates the need to manually enter height for computing BMI. The 350° swivel head allows for line-of-sight viewing and patient privacy.

### Advantages

- ✓ EMR Connectivity**  
 Scale interfaces with a computer, EMR software or other electronic device via Health o meter® Professional Connectivity Solutions
- ✓ 350° Swivel Display**  
 Eye-level rotating display can be read from almost any angle and allows for ease of use and patient privacy
- ✓ Built-in Digital Height Rod**  
 The easy-to-use, state of the art digital height rod provides measurements from 23 5/8" - 88" / 60 cm - 222 cm
- ✓ Aluminum Corrosion-Resistant Column**
- ✓ Value-Added Features**  
 Operating functions (LB/KG Conversion, LB/KG Lock Out, BMI, Zero, Tare, Hold/Release, Auto Zero, Auto Off) included at no additional cost
- ✓ Calculates Body Mass Index (BMI)**  
 Scale functions include the ability to calculate the patient's BMI
- ✓ Patented Design**  
 The 600KL's patented design provides a user-friendly and intuitive interface for measuring patient weight, height and BMI with a convenient rotating display head



### Specifications



US Patent #D667323

- Capacity: 660 lb / 300 kg
- Resolution: 0.2 lb / 0.1 kg
- Height Rod Resolution: 0.1" / 0.1cm
- Platform Size (w x d x h):  
14 1/4" x 14 1/4" x 2 5/8"  
362 mm x 362 mm x 67 mm
- Product Footprint (w x h x d):  
14 3/4" x 24 1/8" x 54 3/8"  
375 mm x 613 mm x 1381 mm
- Product Weight: 28 lb / 13 kg
- Digital Height Rod: 23 5/8" - 88" / 60 cm - 222 cm
- Connectivity: USB, optional Pelstar wireless technology
- 1" / 25 mm LCD display, 35° left / right swivel  
(350° with height rod fully retracted) for ease of viewing
- Functions: LB / KG Conversion, LB / KG Lock Out, BMI, Zero, Tare, Hold / Release, Auto Zero, Auto Off
- 2 wheels for ease of movement
- 100 - 240V adapter (ADPT50) included, 6 C-cell batteries (not included)
- 2 year limited warranty

Model	Description
600KG	Heavy duty eye level digital scale with integral digital height rod, KG only
600KL-BT	600KL with built-in Pelstar wireless technology
600KG-BT	600KG with built-in Pelstar wireless technology
SS-600KL	ScaleSurance 2 year extended warranty for 600KL & 600KG*
SS-600KL-BT	ScaleSurance 2 year extended warranty for 600KL-BT & 600KG-BT*
ADPT50	Adapter included, item # for replacement or spare
PLUGSET50	Plug adapters for use outside North America

Scale also includes connectivity capabilities. To learn about connectivity options, visit [homscales.com/ConnectivitySolutions](http://homscales.com/ConnectivitySolutions)

\* Not available in all countries

600KL\_F&B\_20250513